

Mobile 50-litre clinical waste bin with foot-pedal lid and castors. Epoxy-coated steel for hands-free, hygienic waste collection in wards.

This 50-litre clinical waste bin is fabricated from epoxy-coated steel on a mobile castor base, with a foot-pedal mechanism that raises the lid hands-free to support infection control. The smooth, chip-resistant coating wipes down easily and resists daily chemical cleaning, while the enclosed body contains odour and holds a standard removable liner for ward, theatre and treatment-room use.

SKU: VM-50CB