



The 6-Drawer, 3-Shelf Anaesthetic Trolley is an extended clinical-grade storage trolley built from 1.2mm epoxy powder-coated mild steel with a stainless steel work surface. Combining six organised drawers with three open shelves, it provides significantly more storage capacity than a drawer-only unit ? making it ideal for busy operating theatres and procedure rooms requiring both concealed and open-access equipment storage.

The **6-Drawer, 3-Shelf Anaesthetic Trolley** is an extended-capacity clinical storage trolley designed for operating theatres, anaesthetic rooms, procedure suites, and critical care environments where greater storage flexibility is required. Building on the functionality of a standard 6-drawer anaesthetic trolley, this model adds three open shelves alongside the six drawers ? providing a comprehensive, dual-access storage system for anaesthetic drugs, consumables, monitoring equipment, and larger clinical items that may not fit within standard drawers.

Six Drawers Plus Three Open Shelves

The combination of six drawers and three open shelves gives clinical staff the best of both storage formats. Drawers provide secure, organised concealed storage for medications, syringes, airways, and small consumables, while the open shelves offer rapid, unobstructed access to larger items such as breathing circuits, IV bags, monitoring accessories, and equipment bags ? reducing preparation time and keeping the workspace organised throughout procedures.

Stainless Steel Work Surface

The stainless steel top surface provides a hygienic, non-porous workspace for preparing anaesthetic equipment and medications. It withstands repeated cleaning with hospital-grade disinfectants and corrosive cleaning agents without degrading ? maintaining the hygiene standards essential in surgical environments.

Robust Construction &

SKU: VM-AT6D3S