



The FS793S 360° Swivel Shower Chair (VM-SCS793S) is a swivel shower chair with a 360-degree rotating seat mechanism, designed to assist users transferring from a wheelchair, bath edge or shower chair into position within the shower without standing, lifting or lateral weight-shifting. It has a fixed outer frame measuring 77 × 56 × 38.5cm with non-slip rubber feet, a rotating inner seat, a weight capacity of 100kg, a net weight of 8.8kg and a gross weight of 10.6kg. It is packed 2 pieces per carton. Available in South Africa from Vital Medical.

360° Swivel Shower Chair for Safe Transfer and Seated Bathing

The FS793S is a swivel shower chair with a 360-degree rotating seat mechanism, designed for users who require the additional transfer assistance of a rotating seat during the move from wheelchair, bath edge or shower chair into position within the shower. The swivel mechanism allows the user to sit on the chair at the shower entrance or bath edge and then rotate into the optimal showering position without standing, lifting or lateral weight-shifting, significantly reducing the physical effort and fall risk of the transfer process. For users with severe lower limb weakness, hemiplegic gait patterns or post-arthroplasty precautions, entering the shower and positioning without a standing transfer is a clinically significant functional improvement that supports showering independence or reduces the physical demands on assisting carers.

Stable Base Frame with Rotating Seat

The FS793S separates the functions of stability, provided by the fixed outer frame, and positioning, provided by the rotating inner seat, so that the frame remains stationary while the seat rotates the user smoothly into position. The frame dimensions of 77 × 56 × 38.5cm fit within standard shower cubicle and wet room access dimensions. The 100kg weight capacity and 8.8kg net weight make the FS793S practical for both user and carer handling in daily use. Non-slip rubber feet on the outer frame maintain floor contact during the rotational movement, preventing the chair from sliding during transfers.

South African Transfer-Assist Shower Equipment

The FS793S swivel shower chair is a recommended item in South African occupational therapy and rehabilitation prescriptions for users who cannot perform a standing shower entry or who require the additional support of a rotational transfer mechanism. Available in South Africa from Vital Medical, the FS793S is supplied to South African healthcare facilities, occupational therapists, home care services and individuals requiring a swivel shower chair for safe daily showering with reduced transfer effort. Facilities can request a free quote through the Vital Medical website and receive a personalised response within 24 hours.

TECHNICAL SPECIFICATIONS TABLE

Specification	Detail
Brand	Vital Medical
SKU	VM-SCS793S
Model	FS793S
Type	360° swivel shower chair
Dimensions	77 × 56 × 38.5cm
Seat Rotation	360° swivel mechanism
Weight Capacity	100kg
Rubber Feet	Non-slip
Net Weight	8.8kg
Gross Weight	10.6kg
Packing	2 pieces per carton

SKU: VM-SCS793S