



The Vital Medical Plastic Bedpan (VM-BPP) is a durable clinical-grade plastic bedpan for bedridden patients unable to use a standard toilet. Lightweight, smooth-surfaced for hygiene, and designed for comfortable use in hospital wards, nursing homes, and home care settings.

What Is a Clinical Plastic Bedpan?

The VM-BPP plastic bedpan is a fundamental patient care item for bedridden individuals who are unable to access a toilet ? whether due to acute illness, post-operative recovery, fracture, neurological impairment, or end-of-life care. The bedpan is placed under the patient in the supine or slightly elevated position, allowing urination and defecation in bed without requiring the patient to be transferred to a commode or bathroom.

Clinical and Nursing Application

Bedpans are a routine nursing care item in hospital wards, rehabilitation units, aged care facilities, and home nursing. Proper bedpan technique ? using a protective pad, warming the pan before use, and supporting the patient's back ? minimises discomfort and maintains patient dignity. After use,

the bedpan is removed, covered, and taken to the sluice room for emptying, cleaning, and disinfection before reuse. The smooth plastic surface of the VM-BPP resists staining and supports thorough cleaning with standard healthcare disinfectants.

Material and Advantages Over Stainless Steel

Plastic bedpans offer several advantages over traditional stainless steel: they are lighter weight (easier for nursing staff to handle), they do not conduct heat (warmer against patient skin), and they are less expensive to replace when damaged. The VM-BPP is manufactured from durable clinical-grade plastic that withstands repeated disinfection cycles and normal ward handling. Available alongside the stainless steel bedpan option in the Vital Medical range.

Technical Specifications

Specification		Detail
Specification	Value	
Material	Clinical-grade plastic	
Type	Standard bedpan (supine use)	
Colour	Blue (standard clinical)	
Cleaning	Surface disinfectant compatible	
Application	Bedridden patients ? urination and defecation	
Setting	Hospital ward, aged care, rehabilitation, home nursing	
Complement	Stainless steel bedpan also available	

SKU: VM-BPP