



The FS7931 is a height-adjustable bath bench with handle, 71.5cm length, designed to bridge across a standard bath to allow seated bathing for users who cannot stand or step over the bath edge. Maximum weight capacity 100kg, NW 28kg, GW 33kg. Packed dimensions 32?7?71.5cm. Available from Vital Medical South Africa.

Height-Adjustable Bath Bench with Handle for Safe Seated Bathing

The FS7931 is a height-adjustable bath bench with integrated handle, designed to span across a standard bath and provide a stable seated bathing surface for users who cannot safely stand in the bath or step over the bath edge. The bench allows the user to sit on the rim of the bath at the side, slide across the bench surface into position over the bath, and bathe in a seated position with back support and a secure handle grip ? a configuration that significantly reduces the risk of bath-related falls, which are one of the most common causes of serious injury in older adults and people with mobility impairments. The adjustable legs allow the bench height to be matched to the user's specific bath rim and body height, ensuring a stable, level seating position during bathing.

Safe Bath Transfer for Older Adults and Mobility-Impaired Users

The FS7931 bath bench is particularly suited to older adults, post-surgical patients, and users with lower limb weakness, balance impairments or arthritis who retain some upper body strength and are able to perform a seated lateral transfer. The integrated handle provides a secure grip point for

steadying the transfer and maintaining seated balance during bathing, reducing the reliance on a carer for physical support during the bathing process. The bench surface is designed to be non-slip and comfortable for seated use over the period of a bath.

South African Home Care Bath Safety Equipment

The FS7931 is a practical, effective bath safety solution for South African home users and care services, enabling safe bathing independence for users with moderate mobility impairment without the cost of bath installation modifications such as walk-in shower conversions. Vital Medical supplies the FS7931 to South African home care services, occupational therapists and individuals seeking a safe, affordable bath transfer and seated bathing solution.

Technical Specifications

Specification		Detail
Model	FS7931	
Type	Height-Adjustable Bath Bench with Handle	
Length	71.5cm	
Packed Dimensions	32?7?71.5cm	
Adjustable Legs	Yes	
Weight Capacity	100kg	
Net Weight	28kg	
Gross Weight	33kg	

SKU: VM-BBHA