

The Standard Aneroid Pocket Blood Pressure Meter is an economy-type sphygmomanometer offering reliable 0 to 300 mmHg measurement with a navy blue velcro adult nylon cuff, 2-tube bladder and standard latex bulb and valve – all housed in a vinyl zipper case for easy portability. CE and ISO certified. A cost-effective BP monitoring solution supplied by Vital Medical in South Africa.

### **Reliable Aneroid Blood Pressure Measurement at an Accessible Price Point**

The Standard Aneroid Pocket Blood Pressure Meter is an economy-grade sphygmomanometer designed to deliver accurate blood pressure measurement across the standard clinical range of 0 to 300 mmHg without the cost of premium specialist instruments. Economy-type aneroid gauges use the same fundamental measurement principle as higher-specification models – a pressure-sensitive aneroid capsule deflects a needle across a calibrated dial – making this unit suitable for standard clinical blood pressure assessment in general practice, community health, occupational health and non-critical monitoring environments. CE and ISO certification confirms that the unit meets established quality and calibration standards for medical devices.

### **Velcro Nylon Cuff, 2-Tube Bladder and Vinyl Zipper Case**

The unit includes a navy blue adult-size velcro nylon cuff that fastens quickly and securely without requiring the D-ring winding technique needed for non-velcro cuffs. This makes application faster and reduces the chance of incorrect cuff positioning during routine screening. The 2-tube adult cuff connects to a standard latex inflation bulb and valve, providing controlled pressure inflation and smooth deflation during auscultation. The complete set is stored in a vinyl zipper case measuring 195 × 110 × 75mm and weighing 550 grams – compact and lightweight enough to carry in a bag or jacket pocket for mobile clinical work.

### **Cost-Effective BP Monitoring for South African Health Programmes**

The Standard Aneroid Pocket Blood Pressure Meter is well suited to South African community health screening programmes, school health services and primary care environments where budget constraints require affordable, functional clinical equipment without compromising on calibration

quality. The velcro cuff design speeds up measurement during high-volume screening events, and the compact case ensures the unit is protected between uses. Vital Medical supplies this economy sphygmomanometer to South African healthcare providers, making reliable blood pressure monitoring accessible across a wide range of public and private sector clinical settings.

---

## Technical Specifications

| Specification      | Detail                                  |
|--------------------|-----------------------------------------|
| Type               | Economy Aneroid Pocket Sphygmomanometer |
| Measurement Range  | 0-300 mmHg                              |
| Cuff               | Adult size, navy blue velcro nylon      |
| Cuff Tubes         | 2-tube design                           |
| Inflation          | Standard latex bulb and valve           |
| Case               | Vinyl zipper case                       |
| Packing Dimensions | 195 x 110 x 75mm, 550g                  |
| Certifications     | CE, ISO                                 |

**SKU:** VM-BPMAD