

The Palm Elite Black Aneroid Sphygmomanometer features a 65mm rubber-ringed gauge, non-latex PVC inflation bulb and navy blue adult nylon cuff ? making it a latex-free blood pressure measurement solution for clinicians and patients with latex sensitivities. Calibrated from 0 to 300 mmHg. CE and ISO certified. Available from Vital Medical in South Africa.

Latex-Free Design for Clinicians and Patients with Sensitivities

The Palm Elite Black Aneroid Sphygmomanometer is specifically designed to address the growing need for latex-free clinical equipment in healthcare environments. The inflation bulb is manufactured from PVC rather than latex, eliminating the risk of latex allergen exposure for both the clinician handling the instrument and the patient undergoing blood pressure measurement. Latex allergy is an occupational health concern for healthcare workers with repeated glove and equipment exposure, and it can cause reactions in patients with Type I latex hypersensitivity. The Palm Elite Black allows accurate blood pressure measurement to proceed safely in latex-sensitive environments without compromising measurement quality or clinical workflow.

65mm Rubber-Ringed Gauge and Navy Blue Adult Cuff

The gauge features a 65mm diameter palm-style aneroid face with a protective rubber ring surround that cushions the instrument against drops and knocks, extending gauge service life. The calibration spans 0 to 300 mmHg for standard clinical blood pressure measurement using the Korotkoff auscultation method. The navy blue adult-size nylon cuff is washable and designed to fit standard adult arm circumferences, providing consistent occlusion and reliable readings across a clinical population. The complete unit meets CE and ISO quality and calibration standards.

Versatile Portable Blood Pressure Meter for South African Practice

The Palm Elite Black suits GP practices, hospital wards, occupational health departments and community health programmes in South Africa where latex-free protocols are required or recommended. Its palm-style design is comfortable to hold during auscultation, and the lightweight construction makes it easy to carry between clinical areas. Vital Medical supplies the Palm Elite Black

Aneroid Sphygmomanometer to healthcare providers across South Africa, offering a reliable, latex-safe blood pressure monitoring solution for clinicians committed to patient and staff safety.

Technical Specifications

Specification	Detail
Type	Palm Aneroid Sphygmomanometer (Elite Black)
Dial Diameter	65mm with rubber ring surround
Measurement Range	0-300 mmHg
Cuff	Adult size, navy blue washable nylon
Inflation Bulb	Non-latex PVC
Latex Free	Yes
Certifications	CE, ISO

SKU: VM-BPAPEBK