



The Palm 5 Cuff Aneroid Sphygmomanometer Set includes a precision chrome 65mm diameter dial calibrated from 0 to 300 mmHg, paired with five interchangeable cuffs to accommodate a wide range of patient arm sizes from paediatric to large adult. CE and ISO certified. Ideal for clinics, hospital wards and GP practices requiring multi-patient versatility. Available from Vital Medical.

### **Five-Cuff Versatility for Every Patient Size**

The Palm 5 Cuff Aneroid Sphygmomanometer is designed for clinical settings where accurate blood pressure measurement must be taken across a wide range of patient arm sizes within a single practice or ward environment. The set includes five interchangeable cuffs covering paediatric through to large adult arm circumferences, eliminating the need to maintain separate sphygmomanometers for different patient populations. Correct cuff sizing is a fundamental

requirement for accurate blood pressure measurement ? an oversized or undersized cuff introduces systematic measurement error ? and this multi-cuff set enables clinicians to select the correct cuff for each patient quickly and confidently.

### **65mm Chrome Dial with Precision Calibration**

The unit features a 65mm diameter chrome-plated aneroid gauge calibrated from 0 to 300 mmHg. The palm-style handheld form factor places the gauge directly in the clinician's grip for clear reading while auscultating with a stethoscope in the other hand. The gauge's chrome finish and ABS gauge housing are durable under regular clinical handling, and the calibration is maintained in accordance with CE and ISO quality standards. The palm design is widely used in hospital wards and outpatient departments where clinicians take repeated measurements throughout a shift.

### **Clinical Applications Across South African Healthcare**

This five-cuff set is particularly practical for South African clinics and hospitals that serve diverse patient populations, including children, adults and bariatric patients, without the budget or storage space to maintain multiple single-cuff instruments. GP practices, hospital wards, screening programmes and community health centres will find the multi-cuff system efficient and economical. Vital Medical supplies this sphygmomanometer set to healthcare providers across South Africa, supporting accurate blood pressure monitoring in both public and private healthcare settings.

---

### **Technical Specifications**

| <b>Specification</b> | <b>Detail</b>                              |
|----------------------|--------------------------------------------|
| Type                 | Palm Aneroid Sphygmomanometer ? 5 Cuff Set |
| Dial Diameter        | 65mm chrome-plated                         |
| Measurement Range    | 0?300 mmHg                                 |
| Number of Cuffs      | 5 (paediatric to large adult)              |

| Specification  | Detail                  |
|----------------|-------------------------|
| Gauge Housing  | ABS with chrome finish  |
| Inflation      | Standard bulb and valve |
| Certifications | CE, ISO                 |

**SKU:** VM-BPAM5C