



The Aneroid Palm Single-Handed BP Meter with ABS Gauge is available in blue, featuring a powder-coated ABS aneroid gauge, matching navy blue washable adult nylon cuff, 1-tube design and colour co-ordinated vinyl zipper case. Calibrated 0-300 mmHg. Standard latex bulb and valve. CE and ISO certified. Supplied by Vital Medical in South Africa.

Colour-Coded Palm Aneroid for Single-Handed Clinical Use

The Aneroid Palm Single-Handed Blood Pressure Meter with ABS Gauge is a compact, palm-type sphygmomanometer designed for one-handed blood pressure measurement, available in a distinctive blue colour configuration. The blue powder-coated ABS gauge housing is matched to a navy blue washable adult-size nylon cuff, latex inflation bulb and colour co-ordinated vinyl zipper case, creating a complete, visually unified measurement kit. The palm-type form factor combines the gauge, bulb and inflation valve in a single hand-held unit, enabling the user to inflate, measure and

deflate the cuff using one hand ? a practical advantage for self-monitoring, solo clinical assessments and situations where a clinical assistant is unavailable.

ABS Gauge, 1-Tube Design and Colour-Matched Accessories

The aneroid gauge is housed in a durable ABS (acrylonitrile butadiene styrene) casing, calibrated from 0 to 300 mmHg for the full clinical blood pressure measurement range. The 1-tube cuff design connects to the adult-size cuff and bladder, and the standard latex inflation bulb and valve provide controlled pressure inflation and smooth deflation during the auscultation cycle. All accessories ? gauge, cuff, bulb, case ? are colour co-ordinated in blue, making this unit easy to identify in a multi-unit clinical environment, reducing the risk of equipment mix-up in clinic storerooms and providing an aesthetically distinctive appearance for patient interaction.

Portable Blue ABS Aneroid for South African Clinical and Home Use

The blue ABS palm aneroid weighs 400g in its case (195 ? 110 ? 75mm), making it light and portable for daily clinical use, home visits and community health programmes. CE and ISO certified, the unit meets established quality and calibration standards for medical measurement equipment. Vital Medical supplies this colour-coordinated palm aneroid sphygmomanometer to South African healthcare providers and home users who prefer a visually distinctive, single-handed, portable blood pressure measurement solution.

Technical Specifications

Specification	Detail
Type	Palm Aneroid Sphygmomanometer (Single-Handed, ABS Gauge)
Colour	Blue (powder-coated ABS gauge)
Measurement Range	0?300 mmHg
Cuff	Navy blue washable adult nylon, 1-tube design

Specification	Detail
Inflation	Standard latex bulb and valve
Case	Navy blue vinyl zipper case (colour co-ordinated)
Packing Dimensions	195 ? 110 ? 75mm, 400g
Certifications	CE, ISO

SKU: VM-BPAPABS