



The BP Meter Aneroid Pocket Standard (VM-BPMAPKS) is an economy-grade portable aneroid sphygmomanometer with a calibrated gauge range of 0–300 mmHg, a washable navy blue adult velcro nylon cuff in a 2-tube configuration, and a standard latex bulb and valve. CE and ISO certified, it is supplied complete in a vinyl zipper carry case for convenient storage and transport. The VM-BPMAPKS is suited to general practitioners, nurses, community health workers, pharmacies, first aid kits, and home monitoring where a reliable portable blood pressure monitor is required at an economical price.

Reliable Auscultatory Measurement Across the Clinical Range

The BP Meter Aneroid Pocket Standard (VM-BPMAPKS) is an economy-grade portable aneroid sphygmomanometer designed for accurate, straightforward auscultatory blood pressure measurement across a wide range of clinical and community settings. The aneroid gauge is calibrated from 0 to 300 mmHg, covering the full range required for standard adult blood pressure assessment. The standard latex bulb and valve provide smooth, controlled inflation and deflation for comfortable and repeatable measurements session to session. CE and ISO certification confirms compliance with applicable quality and safety standards for medical measuring instruments.

Cuff Design and Hygienic Reusability

The VM-BPMAPKS is fitted with an adult-size navy blue nylon cuff with velcro fastening, allowing quick application and a secure fit across a range of adult arm circumferences without the need for buckle adjustment. The 2-tube adult cuff and bladder configuration ensures consistent pressure delivery from bulb to cuff during inflation and controlled release during auscultatory measurement. The nylon cuff is washable and reusable, supporting basic infection control practices in clinical and community health settings where multiple patients may be assessed within a single session.

Portability, Storage, and Clinical Applications

The BP Meter Aneroid Pocket Standard is supplied complete in a compact vinyl zipper case that keeps the gauge, cuff, bulb, and valve together and protected between uses. Its lightweight, pocket-sized form factor makes it straightforward to carry in a clinical bag, home visit kit, or emergency pack without adding meaningful bulk or weight. Available in South Africa from Vital Medical, the VM-BPMAPKS is suited to general practitioners, nursing staff, community health workers, pharmacies, occupational health screenings, first aid kits, and home monitoring environments requiring a portable aneroid sphygmomanometer at an accessible price point.

Technical Specifications

Specification	Value
VM SKU	VM-BPMAPKS
Type	Economy aneroid sphygmomanometer
Gauge range	0–300 mmHg
Cuff type	Adult size, navy blue, washable velcro nylon
Cuff configuration	2-tube adult cuff and bladder
Bulb and valve	Standard latex
Carry case	Vinyl zipper case
Certifications	CE; ISO

SKU: VM-BPMAPKS