



This aneroid sphygmomanometer and stethoscope combination set provides clinicians with a complete, portable blood pressure measurement solution. The calibrated aneroid gauge reads from 0 to 300 mmHg and is paired with a navy blue washable adult nylon cuff and matching stethoscope. CE and ISO certified. Supplied by Vital Medical in South Africa.

A Complete Blood Pressure Assessment Kit for Clinical Use

This aneroid sphygmomanometer and stethoscope combination set gives clinicians and healthcare practitioners everything needed for accurate auscultatory blood pressure measurement in a single, ready-to-use package. The aneroid gauge is precisely calibrated across a 0 to 300 mmHg measurement range, providing reliable systolic and diastolic readings through the standard Korotkoff sound method used in clinical practice worldwide. The inclusion of a matched stethoscope eliminates the need to source compatible accessories separately, making this set an efficient and cost-effective first choice for GP practices, outpatient clinics and community health workers requiring portable BP assessment capability.

Navy Blue Nylon Cuff and Quality Construction

The set includes a washable adult-size navy blue nylon cuff designed to fit a standard adult arm circumference comfortably and consistently. The cuff is paired with a standard latex inflation bulb and valve, ensuring controlled inflation and smooth, predictable deflation during auscultation. The overall construction of both the aneroid gauge and stethoscope meets CE and ISO quality standards, providing assurance of measurement accuracy and build reliability across repeated clinical use. The washable cuff supports basic hygiene protocols between patients.

Portable and Practical for South African Healthcare Settings

The sphygmomanometer and stethoscope combination set is well suited to the demands of South African primary care, where clinicians frequently need to move between consultation rooms, ward rounds and community screening programmes. The compact, portable format is easy to carry and store, and the complete kit format reduces equipment procurement complexity for practice managers. Vital Medical supplies this combination set to South African healthcare providers, offering a straightforward and reliable blood pressure monitoring solution backed by local distribution and support.

Technical Specifications

Specification	Detail
Type	Aneroid Sphygmomanometer + Stethoscope Set
Measurement Range	0-300 mmHg
Cuff	Adult size, navy blue washable nylon
Inflation	Standard latex bulb and valve
Accessories Included	Matching stethoscope
Certifications	CE, ISO

SKU: VM-BPAM+STETH