



The Wall Swivel Square Face Aneroid Sphygmomanometer features a 14cm square dial with a 180-degree horizontal swivel manometer, calibrated from 0 to 300 mmHg. The swivel mount allows optimal dial positioning for clinician viewing from multiple angles, and the unit includes a navy blue washable adult nylon cuff and PVC coiled tubing. CE and ISO certified. Available from Vital Medical South Africa.

180-Degree Swivel Manometer for Flexible Viewing in Clinical Spaces

The Wall Swivel Square Face Aneroid Sphygmomanometer is engineered to address a common limitation of fixed wall-mounted BP units: the inability to adjust dial orientation for different clinician positions or room layouts. The 14cm square dial is mounted on a horizontal swivel mechanism that rotates 180 degrees, allowing the clinician to position the gauge face for optimal reading regardless of where they are seated or standing relative to the wall mount. This flexibility is valuable in consulting rooms where the layout requires the clinician and patient to sit in varied orientations, or in multi-use clinical spaces where different practitioners use the same mounted unit.

Large Square Dial with Navy Blue Cuff and PVC Coiled Tubing

The 14cm square aneroid dial is calibrated from 0 to 300 mmHg and provides a high-visibility pressure reading across the standard clinical blood pressure measurement range. The navy blue washable adult-size nylon cuff with D-ring closure ensures correct, repeatable positioning on the patient's arm, and the 2-tube design connects to an adult-size bladder for consistent occlusion. PVC coiled tubing prevents tangling during clinical use and maintains a tidy cuff and gauge setup between measurements. The standard latex inflation bulb and spring-loaded inflation valve provide smooth pressure control during the auscultation cycle.

Practical Wall-Mounted BP Solution for South African Clinical Settings

The swivel mount design is particularly practical for South African consulting rooms and clinic bay setups where clinicians see high patient volumes throughout the day. Having a wall-mounted unit with a flexible viewing angle reduces the need to physically reposition the gauge between patients, improving workflow efficiency. The unit is CE and ISO certified and is supplied by Vital Medical to South African healthcare providers, offering a durable, wall-mounted blood pressure measurement solution suited to both public and private sector clinical environments.

Technical Specifications

Specification	Detail
Type	Wall-Mounted Swivel Aneroid Sphygmomanometer
Dial Shape	Square face
Dial Size	14cm (140mm)
Swivel Range	180 degrees horizontal
Measurement Range	0-300 mmHg
Cuff	Adult size, navy blue washable nylon with D-ring
Cuff Tubes	2-tube design
Tubing	PVC coiled

Specification	Detail
Inflation	Standard latex bulb with spring inflation valve
Packing Dimensions	165 ? 195 ? 175mm, 960g
Certifications	CE, ISO

SKU: VM-BPAWMSQ