



The BP Meter Aneroid Wall/Mobile Square (VM-BPMAWMSQ) is a clinical-grade aneroid sphygmomanometer with a 140mm square-face dial and 90° swivel manometer for flexible, multi-angle blood pressure readings. Left/right adjustable and compatible with wall mounting or a mobile stand, it delivers accurate 0–300 mmHg measurement with a washable navy blue adult nylon cuff and standard latex bulb and valve with spring. CE and ISO certified.

Large Square Dial and 90° Swivel for Clear, Flexible Readings

The BP Meter Aneroid Wall/Mobile Square is a professional clinical aneroid sphygmomanometer designed for accurate, consistent blood pressure measurement in fixed clinical environments. Its defining feature is the 140mm square-face dial — significantly larger than a standard pocket gauge — providing a clear, high-visibility readout that reduces reading error in busy clinical settings. The 90° swivel manometer allows the dial face to be angled toward the clinician regardless of mounting

position or patient orientation, eliminating the need to reposition the unit for each patient and reducing physical strain on clinical staff during high-volume measurement sessions. The left/right adjustable configuration ensures the unit can be orientated to suit the clinician's dominant hand and the consulting room layout.

Precision Measurement and Hygienic Cuff Design

The aneroid gauge is calibrated from 0–300 mmHg, covering the full diagnostic range required for routine adult blood pressure assessment. The standard latex bulb and inflation valve with spring deliver smooth, controlled inflation and a consistent deflation rate for a comfortable patient experience and accurate Korotkoff sound identification. The adult-size navy blue washable nylon cuff is reusable and can be laundered between patients, supporting standard infection-control practice in clinical environments where multiple patients are assessed using the same equipment.

Mounting Versatility and Clinical Applications

The Wall/Mobile Square is designed for fixed-station clinical deployment, mounting to a wall bracket for permanent consulting room installation or fitting to a mobile stand for ward-round or multi-room use. This versatility makes it suited to hospital outpatient departments, ward nursing stations, general practice consulting rooms, occupational health clinics, and any clinical environment where a high-visibility, robustly mounted sphygmomanometer is preferred over a handheld or portable unit. CE and ISO certified. Available in South Africa from Vital Medical.

TECHNICAL SPECIFICATIONS TABLE

Specification	Detail
SKU	VM-BPMAWMSQ
Type	Professional aneroid sphygmomanometer
Dial	140mm square face

Specification	Detail
Manometer	90° swivel
Adjustment	Left/right adjustable
Gauge Range	0-300 mmHg
Cuff	Adult size, navy blue washable nylon
Inflation	Standard latex bulb and valve with spring
Mounting	Wall / mobile stand
Certifications	CE, ISO

SKU: VM-BPMAWMSQ