



The Wall/Desk Square Face Aneroid Sphygmomanometer features an oversized 145mm square dial calibrated from 0 to 300 mmHg, designed for high-visibility blood pressure monitoring in consulting rooms, hospital wards and reception areas. The latex-free system includes a navy blue adult nylon cuff and standard inflation bulb. CE and ISO certified. Supplied by Vital Medical in South Africa.

High-Visibility Square Dial for Consulting Room Mounting

The Wall/Desk Square Face Aneroid Sphygmomanometer is designed for fixed mounting in consulting rooms, wards and clinical reception areas where a large, easy-to-read display improves measurement speed and accuracy. The oversized 145mm square dial provides a clear view of the pressure needle from a greater distance than standard pocket models, reducing measurement error caused by parallax or poor lighting in busy clinical environments. The gauge is calibrated from 0 to 300 mmHg and is designed for wall mounting or desk positioning, allowing the clinician to take blood pressure readings while maintaining a comfortable seated position without needing to hold the gauge.

Latex-Free Construction with Washable Adult Cuff

The unit uses a latex-free inflation system, making it suitable for use in environments that require or recommend avoidance of latex allergens for patient and staff safety. The navy blue adult-size nylon cuff is washable, supporting hygiene protocols between patient consultations. The standard inflation bulb and valve provide smooth, controlled inflation and deflation during the auscultation cycle. The 145mm square gauge housing is robustly constructed for long-term fixed installation in a clinical setting, and the overall unit meets CE and ISO quality certification requirements for medical measurement equipment.

Fixed-Station BP Measurement for South African Clinics and Wards

This wall/desk sphygmomanometer is ideal for South African GP practices, occupational health units and hospital outpatient departments that conduct high volumes of blood pressure assessments from a fixed clinical station. The large square dial is particularly useful in lower-light environments or where clinicians take multiple readings quickly during busy clinic sessions. Vital Medical supplies this unit to South African healthcare providers, offering a durable, high-visibility blood pressure monitoring solution for fixed-station clinical use across primary care, specialist practice and hospital environments.

Technical Specifications

Specification	Detail
Type	Wall / Desk Aneroid Sphygmomanometer
Dial Shape	Square face
Dial Size	145mm
Measurement Range	0-300 mmHg
Cuff	Adult size, navy blue washable nylon
Latex Free	Yes
Mounting	Wall or desk

Specification	Detail
Inflation	Standard bulb and valve
Certifications	CE, ISO

SKU: VM-BPAWD2004