

The Pocket Deluxe Aneroid Sphygmomanometer features a shock-proof chrome-finished gauge calibrated from 0 to 300 mmHg, with a washable cotton cuff for comfortable, hygienic blood pressure measurement in clinical and community settings. Its compact pocket design makes it ideal for GPs, nurses and community health workers requiring a reliable portable BP meter. CE and ISO certified. Supplied by Vital Medical.

Shock-Proof Chrome Gauge for Portable Clinical Use

The Pocket Deluxe Aneroid Sphygmomanometer is engineered for clinicians who require a reliable, durable pocket-format blood pressure meter that can withstand the physical demands of daily use in mobile and clinical settings. The gauge features a shock-proof chrome-finished construction calibrated from 0 to 300 mmHg, providing accurate systolic and diastolic readings via the Korotkoff auscultation method. The shock-proof design protects internal calibration from the drops and impacts that inevitably occur in busy ward environments, outpatient clinics and community health programmes, maintaining measurement accuracy over an extended service life.

Washable Cotton Cuff for Patient Comfort and Hygiene

Unlike standard nylon cuffs, the Pocket Deluxe uses a washable cotton cuff that is softer against the skin and comfortable for patients with sensitive skin or latex sensitivities. The cotton cuff can be removed and laundered between patients, supporting basic infection control and hygiene protocols in multi-patient environments. The cuff is designed to seat correctly on a standard adult arm, ensuring consistent occlusion and reliable blood pressure readings when inflated via the standard latex bulb and valve system.

Practical for South African Community and Primary Health Care

The compact pocket format of this sphygmomanometer suits the working reality of South African community health workers, school nurses, occupational health practitioners and primary care nurses who conduct blood pressure screening across multiple sites in a single day. The instrument is light, easy to carry in a pocket or bag, and robust enough for daily clinical use in resource-limited settings.

CE and ISO certified, this unit is supplied to South African healthcare providers by Vital Medical, delivering dependable blood pressure measurement capability across primary, community and occupational health settings.

Technical Specifications

Specification	Detail
Type	Pocket Aneroid Sphygmomanometer (Deluxe)
Measurement Range	0-300 mmHg
Gauge	Shock-proof chrome-finished aneroid
Cuff	Washable cotton, adult size
Inflation	Standard latex bulb and valve
Certifications	CE, ISO

SKU: VM-BPAD