



The MS7301 is a digital adult floor scale with a 250kg capacity at 100g increments, automatic step-on activation and auto-off, wide 410 x 360mm PVC platform for bariatric patients, self-adjusting feet, non-slip surface, and a 1.0" LCD display. CE/ISO certified, powered by 4 x AAA batteries, and designed for reliable professional-grade weighing in clinical, occupational health, and high-traffic screening environments.

The **MS7301 Digital Adult Floor Scale** is a professional-grade floor scale combining a high-capacity 250kg measurement range with the practical features required in clinical and occupational health environments – including automatic step-on activation, a wide bariatric-accommodating platform, and a low-profile design for easy access.

### **250kg Capacity x Auto Step-On Activation**

The MS7301 activates automatically when a patient steps onto the platform and shuts off automatically after use – eliminating the need to press any buttons and streamlining high-volume weighing workflows. The 0-250kg capacity at 100g increments covers the full adult weight range including higher-weight patients.

### **Wide Platform – Bariatric Accommodation**

The 410 x 360mm PVC surface platform is wider than standard consumer scales – providing the stability and standing space required for bariatric patients and for patients with balance or mobility limitations. The non-slip PVC surface texture reduces slip risk during measurement.

### **Low Profile – Easy Access**

The low-profile platform design minimises the step height required to mount the scale – an important accessibility feature for patients with limited mobility, lower limb weakness, or balance impairment, and for any environment where rapid patient throughput is required.

### **Self-Adjusting Feet – Anti-Tip Stability**

Self-adjusting levelling feet automatically compensate for minor floor irregularities, keeping the scale stable and level on uneven surfaces and preventing tipping during use.

### **CE/ISO Certified**

The MS7301 carries CE and ISO quality standard certification, confirming compliance with the relevant health, safety, and performance standards for medical measuring equipment.

**Suitable For:** GP practices, hospital outpatient departments, bariatric clinics, occupational health

assessments, pharmacy wellness stations, and high-traffic health screening environments.

## Technical Specifications

| Specification                    | Detail                                |
|----------------------------------|---------------------------------------|
| Product Code                     | SCA7301                               |
| Model                            | MS7301                                |
| Capacity                         | 0 ? 250kg                             |
| Increment                        | 100g                                  |
| Display                          | 1.0" LCD                              |
| Platform Size                    | 410 ? 360mm (PVC surface)             |
| Overall Dimensions (with carton) | 460 ? 400 ? 75mm                      |
| Scale Weight                     | 3.5kg                                 |
| Activation                       | Automatic step-on / auto-off          |
| Feet                             | Self-adjusting (anti-tip)             |
| Surface                          | Non-slip texture                      |
| Power                            | 4 ? 1.5V AAA batteries (not included) |
| Quality Standard                 | CE / ISO                              |

**SKU:** VM-MS7301