



The Charder Cupid 1 Digital Baby Scale (VM-SBC1), model SCAC1, is a digital baby scale with a removable tray designed for accurate infant weight measurement in home, hospital and infant care centre settings. It has a 0-20kg capacity, with 10g increments across the 0-10kg range, an overall footprint of 320 × 525 × 160mm with a 320 × 520mm platform, a 1.0" LCD display, and function keys for ON/OFF, KG/LB switching, ZERO, HOLD and TARE. It weighs 2.2 kg and runs on a 9V alkaline battery or 230V AC adaptor. It is CE certified. Available in South Africa from Vital Medical.

Digital Baby Scale for Home, Hospital and Clinic Use

The Charder Cupid 1 (VM-SBC1, SCAC1) is a digital baby scale with a removable tray, designed for accurate infant weight measurement in home, hospital and infant care centre settings. The 0-20kg capacity covers the full infant and toddler weight range, with 10g increments providing clinically appropriate precision across the 0-10kg range, where small weight changes are among the most significant clinical indicators of infant feeding adequacy and growth progress. The pan-style tray provides a safe, contained weighing surface for infants who cannot yet support themselves, and the tray is removable for older babies and toddlers who can stand on the base platform, making the Cupid 1 a versatile scale that remains appropriate across the full growth trajectory from newborn to toddler.

Practical Dimensions and Controls

The Cupid 1's 320 × 525mm platform provides ample space for positioning an infant safely within the tray. The 1.0" LCD display is clearly readable by a parent or clinician positioned beside the scale. Function keys include ON/OFF, KG/LB switching, ZERO, HOLD and TARE. The HOLD function is the most practically important for infant weighing, as it locks the displayed measurement when the infant is moving, preventing the need for repeated weighing events to capture a stable reading. At 2.2 kg, the Cupid 1 is lightweight and portable for use in both fixed clinical settings and home visits.

South African Infant and Toddler Weighing Scale Supply

The Cupid 1 is CE certified and powered by a 9V alkaline battery or a 230V AC adaptor, providing flexible power options for both clinical and home use. Available in South Africa from Vital Medical, the Cupid 1 is supplied to South African maternity hospitals, paediatric clinics, community health centres and home users requiring a reliable, accurate baby scale for routine infant weight monitoring. Facilities can request a free quote through the Vital Medical website and receive a personalised response within 24 hours.

TECHNICAL SPECIFICATIONS TABLE

Specification	Detail
Brand	Charder
SKU	VM-SBC1
Model	Cupid 1 (CODE: SCAC1)
Type	Digital baby scale with removable tray
Capacity	0-20kg (0-10kg: 10g increments)
Overall Dimensions	320 × 525 × 160mm
Platform Size	320 × 520mm
Display	1.0" LCD
Power Supply	9V alkaline battery / 230V AC adaptor
Unit Weight	2.2kg
Removable Tray	Yes
Function Keys	ON/OFF, KG/LB, ZERO, HOLD, TARE
Quality Standard	CE

SKU: VM-SBC1