

The FS691 is a powder-coated steel commode chair with detachable armrests, detachable footrests, locking rear castors and a plastic commode seat. Seat width 44cm, seat height 38cm, overall height 88cm, 100kg capacity, NW 12kg. For bedside and bathroom toileting use. Available from Vital Medical South Africa.

### **Steel Commode Chair with Detachable Arms, Footrests and Locking Castors**

The FS691 is a powder-coated steel commode chair that combines the functionality of a wheeled transport chair with the clinical utility of a bedside or bathroom commode. The chair is designed for users who are unable to walk to the toilet or who require regular toileting assistance in a bedside or ward environment. The detachable armrests and footrests allow user transfer from bed, chair or wheelchair to the commode without needing to navigate fixed side components – a critical safety feature for reducing transfer-related injury risk and enabling carers to assist safely from either side. The solid rear castors with locking mechanism allow the chair to be wheeled into position and then locked firmly before the user transfers or uses the commode seat, preventing chair movement during use.

### **Practical Commode Dimensions for Ward and Home Use**

The FS691's 44cm seat width, 38cm seat height and 44cm seat depth provide a comfortable, stable seated position for adult users. The 51cm height from seat to floor is compatible with standard bed and transfer heights. The 55cm overall width allows passage through standard doorways and ward bays. The 12.5cm diameter castors provide smooth rolling on hard floors with sufficient castor size to manage minor surface irregularities. The powder-coated steel frame is durable, easy to clean and resistant to the repeated disinfection required in clinical environments.

### **South African Clinical and Home Commode Supply**

The FS691 steel commode chair is a standard item for South African hospital wards, rehabilitation units, frail care facilities and home care services where bedside toileting capability is required for patients and residents with limited mobility. Vital Medical supplies the FS691 to South African

healthcare facilities and home care providers requiring a reliable, functional commode chair for daily clinical or care use.

---

### Technical Specifications

| Specification          | Detail                                              |
|------------------------|-----------------------------------------------------|
| Model                  | FS691                                               |
| Type                   | Steel Commode Chair with Detachable Arm & Footrests |
| Frame                  | Powder-coated steel                                 |
| Seat Width             | 44cm                                                |
| Seat Height            | 38cm                                                |
| Seat Depth             | 44cm                                                |
| Height (seat to floor) | 51cm                                                |
| Overall Length         | 48cm                                                |
| Overall Width          | 55cm                                                |
| Overall Height         | 88cm                                                |
| Castor Diameter        | 12.5cm                                              |
| Rear Wheel Diameter    | 12.5cm                                              |
| Weight Capacity        | 100kg                                               |
| Net Weight             | 12kg                                                |
| Packing (L?W?H)        | 54?53?53cm                                          |

**SKU:** VM-CCFS691