

Cot Bed Child Drop-down Cotsides Incl 75mm Swivel & Brake Castors.
1370x760 Wire mesh Mattress base



Child cot bed with drop-down cotsides, wire-mesh base and braked castors. Safe, mobile paediatric bed, 1370 x 760mm.

This paediatric cot bed features drop-down cotsides for safe, easy patient access over a 1370 x 760mm wire-mesh mattress base that promotes ventilation. Mounted on 75mm swivel castors with brakes, it moves smoothly and locks securely in place. The durable, easy-clean frame is designed for children's wards, paediatric units and hospital nurseries.

SKU: VM-CBCDPD