

The Vital Medical Delivery Kit Pack 10pcs (VM-DK10P) is a pre-assembled sterile obstetric set containing all instruments and consumables required for a controlled vaginal delivery. Includes cord clamps, obstetric scissors, forceps, delivery drape, swabs, and gloves.

What Is the Delivery Kit Pack?

The VM-DK10P is a 10-piece sterile obstetric delivery kit designed for midwives, general practitioners, and paramedics performing controlled vaginal deliveries in hospital, clinic, or community settings. It consolidates the essential instruments and sterile consumables for a safe delivery into one ready-to-use pack, reducing preparation time and minimising the risk of missing a critical item in an urgent delivery situation.

Contents and Delivery Procedure

The kit contains the complete instrument set for neonatal cord management and perineal care: umbilical cord clamps (3), cord cutting scissors, episiotomy/delivery scissors, Kocher haemostatic forceps, sterile delivery drape, maternity swabs, and disposable examination gloves. Following delivery of the infant, the cord is double-clamped, divided, and the perineum assessed for lacerations or episiotomy repair requirements. The 10-piece set covers all routine steps of the third stage through to initial neonatal assessment.

Application Across Settings

The VM-DK10P is suitable for hospital labour wards, midwife-operated obstetric units, community health clinics, and paramedic ALS bags for emergency deliveries. The compact format of a pre-assembled kit is particularly valuable in unplanned or out-of-hospital deliveries where time and resources are limited. Complement with a resuscitator, suction equipment, and warm blankets for complete neonatal care.

Technical Specifications

Specification		Detail
Specification	Value	
Set Contents	10 pieces	
Includes	Cord clamps (3), cord scissors, delivery scissors, Kocher forceps, drape, maternity swabs, gloves	
Application	Controlled vaginal delivery, cord management	
Setting	Hospital, clinic, community, paramedic	
Sterility	Pre-assembled sterile pack	
Reuse	Stainless instruments autoclavable; consumables single-use	

SKU: VM-DK10P