



The Deluxe Dual Head Stethoscope (STJDD/STHDDH) features a silver dual head, 47mm adult head diameter, 54cm tube length, standard black Y-tubing, CE and ISO approved. Single unit weighs 110g. Available in black, blue, green, red and purple. Available from Vital Medical South Africa.

What Is a Dual Head Stethoscope and Why Does It Matter?

A dual head stethoscope has two auscultation surfaces on one chest piece: a diaphragm and a bell. The Deluxe Dual Head Stethoscope (STJDD/STHDDH, VM-STJDD) has a silver dual head with a 47mm adult head diameter, which gives standard acoustic coupling for adult chest wall auscultation. The diaphragm picks up high-pitched sounds such as breath sounds and normal heart sounds. The bell picks up low-frequency sounds such as third and fourth heart sounds and the low-pitched rumble of mitral stenosis. Having both surfaces lets the clinician cover the full range of sounds in a systematic physical examination.

Lightweight Design and Tubing for Daily Clinical Use

At 110g, the STJDD is a lightweight instrument, practical for long clinic sessions, ward rounds and home visits. The standard Y-tubing is 54cm long, which allows comfortable positioning during auscultation without the acoustic degradation associated with excessively long tubes. The black Y-tubing carries sound through the binaural tubing assembly to the earpieces. The single unit is packed in a compact 190 × 100 × 20mm box, so it is easy to store and carry as a standard item of clinical equipment. CE and ISO approvals confirm compliance with international quality standards for clinical diagnostic instruments.

Supply for South African Doctors, Nurses and Students

The STJDD is a cost-effective, widely used stethoscope for South African general practitioners, doctors, nurses and students who need a reliable dual head instrument for daily examinations. Available in South Africa from Vital Medical, it comes in black, blue, green, red and purple, so clinicians in shared clinical environments can identify their own stethoscope. Healthcare facilities can request a free quote through the Vital Medical website and receive a personalised response within 24 hours.

TECHNICAL SPECIFICATIONS TABLE

Specification	Detail
SKU	VM-STJDD
Model	STJDD / STHDDH
Type	Deluxe dual head stethoscope, adult
Head	Silver dual head
Head Diameter (Adult)	47mm
Tube Length	54cm
Y-Tubing	Standard (black; other colours available)
Weight	110g
Available Colours	Black, blue, green, red, purple
Packing (single unit)	190 × 100 × 20mm
Quality Approvals	CE, ISO

SKU: VM-STJDD