



The DW1060 is a complete wall-mounted diagnostic station comprising a 3.5V coaxial ophthalmoscope, 3.5V fibre otoscope, 485 wall-mount transformer with handles, specula dispenser, aneroid sphygmomanometer, and ear specula in two sizes. Designed for permanent installation in consulting rooms, it consolidates the core diagnostic instruments for eye, ear, and blood pressure examination in a single, mains-powered wall unit.

Complete Diagnostic Station for the Consulting Room

The DW1060 integrates three fundamental clinical tools – ophthalmoscope, otoscope, and sphygmomanometer – into a single wall-mounted system, eliminating the need for separate handheld instruments with disposable batteries. The 485 wall-mount transformer provides continuous mains power to both instruments simultaneously, with dedicated handles that park cleanly on the unit between uses. This eliminates the common clinical frustration of depleted batteries mid-examination and ensures consistent, reliable illumination during every assessment.

3.5V Coaxial Ophthalmoscope and Fibre Optic Otoscope

The 3.5V coaxial ophthalmoscope delivers a bright, focused beam for fundoscopic examination of the retina, optic disc, and vasculature. Coaxial design aligns the light source with the viewing axis, reducing corneal reflection and improving visualisation of posterior segment structures. The 3.5V fibre optic otoscope provides a cool, bright white light via fibre optic transmission – minimising heat at the speculum tip during prolonged ear canal and tympanic membrane examination. Together, these instruments cover the full scope of routine head examination in general practice, paediatrics, emergency, and ENT settings.

Specula Dispenser, Dual-Size Ear Specula, and Sphygmomanometer

The included specula dispenser holds disposable ear specula in both clinical sizes – 2.7mm (paediatric) and 4.2mm (adult) – with 35 pieces of each supplied. Wall-mounting the dispenser at the point of use reduces bench clutter and streamlines examination workflow. The aneroid sphygmomanometer completes the station, enabling blood pressure measurement without a separate device. The wooden mounting plate (middle size) accommodates standard wall installation and provides a professional, integrated appearance suited to any consultation or treatment room.

Technical Specifications:

Specification		Value
Model	DW1060	
Type	Wall-mounted diagnostic set	
Ophthalmoscope	3.5V coaxial	
Otoscope	3.5V fibre optic	

Specification	Value
Wall Transformer	485 series with instrument handles
Power	Mains-powered (wall transformer)
Specula Dispenser	Included
Ear Specula Sizes	2.7mm (paediatric) and 4.2mm (adult)
Ear Specula Quantity	35 pcs per size (70 pcs total)
Sphygmomanometer	Aneroid, included
Mounting Plate	Wooden, middle size
Included Accessories	Power cable with plug, wooden mounting plate, specula dispenser, 70 ear specula (35 ? 2.7mm, 35 ? 4.2mm)
Indications	Fundoscopy, otoscopy, tympanic membrane assessment, blood pressure measurement

SKU: VM-DW1060WM