



The Universal Deluxe Diagnostic Set (DIAGSUDBK) includes stainless steel handle, otoscope, ophthalmoscope, laryngeal mirrors 4&5, tongue blade holder, metal tongue depressor, nasal speculum, ear tip connector, ear tips, black leather case. 240?195?35mm, 880g. Requires 2?1.5V C-cell batteries (not supplied). CE and ISO certified. Available from Vital Medical South Africa.

Universal Deluxe Diagnostic Set with Ear Tips and Black Leather Case

The Vital Medical Universal Deluxe Diagnostic Set (code DIAGSUDBK) is the premium version of the universal portable diagnostic set, adding ear tip connector, ear tips and a black leather carry case to the complete instrument complement of the standard universal set ? providing upgraded presentation and additional accessory functionality in a professional black leather case suitable for consultants, specialists and senior clinical staff who require an instrument set that reflects the quality of their clinical practice. The instrument complement ? stainless steel handle, otoscope,

ophthalmoscope, laryngeal mirrors sizes 4 and 5, tongue blade holder, metal tongue blade depressor and nasal speculum ? provides the complete clinical diagnostic toolkit for ENT, ophthalmological and oropharyngeal examination in portable format.

Premium Leather Case Presentation and Ear Tip Accessories

The black leather carry case (240?195?35mm) provides enhanced instrument protection and professional presentation compared to standard hard PVC cases. The leather case format is appropriate for consultant and specialist practice settings where instrument presentation reflects professional standing. The ear tip connector and ear tips provide the ability to use the otoscope with the clinician's own ear tips, allowing personalised ergonomic preference and improved examination positioning, particularly for clinicians who perform high volumes of ear examinations and prefer the tactile control of an ear tip over a standard speculum in the otoscope head.

South African Specialist and Senior Clinical Diagnostic Instrumentation

The Universal Deluxe Set's 880g weight and compact 240?195?35mm case dimensions make it a practical portable set for senior clinicians, specialists, registrars and medical officers who carry their diagnostic set between consulting rooms, wards and clinics. The stainless steel instruments resist corrosion and are compatible with standard sterilisation procedures. CE and ISO certification confirms compliance with international quality standards. Vital Medical supplies the Universal Deluxe Diagnostic Set to South African consultants, specialists, medical officers, nursing sisters and clinical education programmes.

Technical Specifications

Specification		Detail
Product Code	DIAGSUDBK	
Type	Universal Deluxe Diagnostic Set	

Specification	Detail
Contents	Stainless steel handle; Otoscope; Ophthalmoscope; Laryngeal mirrors sizes 4 & 5; Tongue blade holder; Metal tongue blade depressor; Nasal speculum; Ear tip connector; Ear tips
Case	Black leather carry case
Case Dimensions	240?195?35mm
Unit Weight	880g
Carton	20 units per carton, 400?470?230mm
Power	2?1.5V C-cell batteries (not supplied)
Quality Standard	CE, ISO

SKU: VM-DIAGUNB