



The Diagnostic Set Universal LED ? Blue is a complete handheld examination kit comprising an LED otoscope, LED ophthalmoscope, and universal handle, supplied in a white carry case. The bright LED light source delivers consistent, cool illumination for ear, nose, throat, and eye examinations in clinical and primary care settings.

Complete LED Illumination for ENT and Ophthalmic Examinations

The Diagnostic Set Universal LED ? Blue combines two essential diagnostic instruments in a single coordinated kit. The LED otoscope head enables clear visualisation of the external ear canal and tympanic membrane, while the LED ophthalmoscope head supports anterior and posterior eye assessment. Both instruments share a universal handle, reducing the cost and complexity of maintaining separate devices. The blue instrument finish housed in a white carry case ensures easy identification and a professional presentation in any clinical environment.

Advantages of LED Technology in Diagnostic Instruments

LED light sources have largely replaced halogen bulbs in modern diagnostic sets due to their superior colour rendering, energy efficiency, and extended lifespan. The white-balanced LED beam produced by this set accurately represents tissue colour, supporting more reliable clinical assessment compared to the yellowed output of conventional halogen instruments. LED elements are also shock-resistant, making them well suited to the rigours of daily examination use across general practice, emergency, paediatric, and occupational health settings.

Suitable for General Practice, Clinics, and Medical Training

This diagnostic set is designed for use across a broad range of healthcare settings where both ENT and ophthalmic assessment are required. The universal handle system allows rapid interchangeability between the otoscope and ophthalmoscope heads, streamlining examination workflows. Disposable ear specula in assorted sizes are compatible with the otoscope head, supporting infection control protocols. The carry case keeps all components organised and protected between uses, making the set equally appropriate for fixed clinical rooms and mobile or visiting practitioner use.

Technical Specs Table:

Specification	Value
Set Includes	LED otoscope head, LED ophthalmoscope head, universal handle
Light Source	LED
Instrument Colour	Blue
Case Colour	White
Specula Compatibility	Standard disposable ear specula (assorted sizes)
Handle Type	Universal (shared)
Intended Use	ENT and ophthalmic examination
Setting	General practice, clinic, medical training

SKU: VM-DIAGSUBL