



Height-adjustable 3-leg drip stand in epoxy-coated steel on castors. Compact, mobile IV infusion pole for tighter spaces.

This three-leg drip stand features an epoxy-coated steel column on a compact three-leg castor base, with smooth height adjustment and hooks for IV bags and fluids. The narrower footprint manoeuvres easily in confined spaces, while the chip-resistant finish wipes clean. A practical, economical infusion pole for wards, consulting rooms and clinics.

Technical Specifications:

Epoxy powder coated,
3x 75mm swivel castors,
adjustable T- rod,
adjustable height (1300mm-2000mm).

SKU: VM-DS3LEP