



Heavy-duty 5-leg stainless steel drip stand on castors. Maximum stability for multiple IV bags and pump-mounted infusions.

The **5-Leg Mobile Infusion Pole Drip Stand** is a premium, brushed stainless steel IV infusion stand designed for use in hospital wards, operating theatres, treatment rooms, and procedure areas where safe, reliable delivery of intravenous fluids is required. Built to balance clinical-grade stability with ease of patient mobility, it is an essential piece of equipment for any ward or clinical space administering infusion therapy.

Heavy 5-Legged Base for Anti-Topple Stability

The wide, 5-legged heavy base is specifically engineered to prevent accidental toppling – a critical safety consideration in busy ward environments where the stand may be knocked by passing staff, equipment, or mobile patients. The five-point footprint distributes weight evenly and provides a low centre of gravity, keeping the stand upright and stable even when IV bags are loaded and the stand is in motion.

One-Hand Adjustable Height with Safety Twist Lever

The height-adjustable pole is operated via a smooth one-hand twist lever mechanism, allowing clinical staff to quickly set the IV bag height to the correct level for each patient without requiring tools or assistance. The lever is specifically designed so it cannot fall onto the user's thumb – a purposeful safety feature that protects staff during repeated height adjustments throughout a shift.

Patient-Friendly Mobility

The slim, streamlined profile of the infusion pole is sized so that patients can grip and move the stand themselves when ambulating around the ward, supporting patient independence and mobility during infusion therapy. This reduces the demand on nursing staff and supports faster patient recovery through encouraged movement.

Brushed Stainless Steel Finish

The brushed stainless steel frame is corrosion-resistant, easy to wipe clean, and compatible with hospital-grade disinfectants ? maintaining a hygienic surface through the repeated cleaning cycles required in clinical environments. The finish is durable, professional in appearance, and well suited to both clinical and private ward settings.

Suitable For:

All hospital wards, operating theatres, treatment rooms, day clinics, procedure rooms, recovery areas, and outpatient infusion suites.

Technical Specifications

Specification	Detail
Frame Material	Brushed stainless steel
Base	5-legged heavy anti-topple base
Height Adjustment	One-hand twist lever (thumb-safe design)
Mobility	Patient self-manoeuvrable
Suitable For	Wards, theatres, treatment rooms, operating rooms
<i>Number of Hooks</i>	2 Hooks

SKU: VM-DS5LSSR