



Replacement thermal recording paper for the CONTEC ECG1200G 12-channel ECG machine. 210mm wide 30m long roll, high-speed thermal paper compatible with the ECG1200G thermal printing

system. Essential consumable for continuous ECG recording operation. Available from Vital Medical in South Africa.

ECG Recording Paper for CONTEC ECG1200G

The ECG1200G thermal recording paper is the dedicated consumable for the CONTEC ECG1200G 12-channel electrocardiograph, providing the recording medium for all three ECG output modes, auto record, rhythm record and manual record. The 210mm wide roll format accommodates the full 12-lead waveform print layout at the ECG1200G's standard print width, ensuring all leads are reproduced at full scale without trimming or compression. Each 30-metre roll provides extended continuous operation between paper changes, supporting high-volume clinical sessions in busy hospital outpatient departments, ward-based ECG services and cardiac units where recording frequency is high.

Thermal Paper Quality and ECG Waveform Fidelity

Thermal ECG paper quality directly affects the legibility and archival durability of printed ECG records. The ECG1200G paper is a high-speed thermal paper with the surface characteristics required by the ECG1200G's thermal print head to produce sharp, well-defined waveform traces at all supported paper speeds from 5mm/s to 50mm/s and at all sensitivity settings. Clear, consistent waveform reproduction is essential for accurate clinical interpretation, particularly for fine waveform details such as P wave morphology, ST segment displacement and T wave inversion that can be diagnostically important. Properly matched thermal paper also extends the operational life of the thermal print head by maintaining the correct print head-to-paper contact characteristics.

Consumable Supply for South African Cardiac Diagnostic Units

Maintaining an adequate supply of ECG recording paper is essential for uninterrupted cardiac diagnostic services in South African hospitals, outpatient clinics, cardiology units and GP practices operating the CONTEC ECG1200G. Vital Medical supplies the ECG1200G replacement paper as a standalone consumable, enabling South African healthcare facilities to order paper alongside or

independently of the ECG machine itself. Ordering through the equipment supplier ensures paper specification compatibility and consistent waveform quality across the ECG1200G's operational life.

Technical Specifications

Specification	Detail
Type	Thermal ECG Recording Paper
Compatible Device	CONTEC ECG1200G
Width	210mm
Roll Length	30m
Paper Type	High-speed thermal paper
Format	Roll
Application	Auto record, rhythm record, manual record

SKU: VM-ECG1200P