



The Vital Medical Episiotomy Pack 6pcs (VM-EP6) is a pre-assembled stainless steel instrument set for performing episiotomy and perineal repair during obstetric deliveries. Contains Rampley forceps, Braun-Stadler episiotomy scissors, Gillies tissue forceps, Mayo-Hegar needle holder, dental syringe, and operating scissors.

What Is the Episiotomy Pack?

The VM-EP6 is a 6-piece stainless steel obstetric instrument set containing the essential instruments for performing an episiotomy (a controlled perineal incision to facilitate delivery) and subsequent primary repair of the episiotomy wound. Having all required instruments pre-assembled in one sterile set reduces setup time and eliminates the risk of missing a critical instrument during a delivery that may require urgent episiotomy.

Contents and Clinical Use

The set contains: Rampley sponge holding forceps 25cm (for swabbing and maintaining a sterile

field), Braun-Stadler episiotomy scissors 14.5cm (specifically designed for the mediolateral or midline episiotomy incision), Gillies tissue forceps 15cm 1?2 teeth (for wound edge handling during repair), Mayo-Hegar needle holder 18cm (for suturing the episiotomy wound), dental syringe (for local anaesthetic infiltration before the procedure), and operating scissors sh/bl straight 16.5cm (for suture cutting during repair). Together these instruments cover the complete episiotomy procedure from anaesthesia through to wound closure.

Quality and Sterilisation

All instruments in the VM-EP6 are manufactured from surgical-grade stainless steel and are autoclavable at 134?C (prevacuum) for reuse between patients. The set should be processed through CSSD after each use, inspected for instrument integrity, and re-assembled in sterile packaging. The VM-EP6 is suitable for hospital labour wards, obstetric units, and midwife-led birthing centres.

Technical Specifications

Specification	Detail
Specification	Value
Set Contents	6 instruments
Includes	Ramsey Sponge Holding Fcps 25cm, Braun-Stadler Episiotomy Scissors 14.5cm, Gillies Tissue Fcps 15cm 1?2T, Mayo-Hegar NH 18cm, Dental Syringe, Op-Scissors sh/bl str 16.5cm
Material	Surgical-grade stainless steel
Sterilisation	Autoclavable (134?C steam)
Application	Episiotomy and perineal repair
Setting	Labour ward, obstetric unit, midwife-led unit

SKU: VM-EP6