



Ratchet-type examination couch with a ratchet type adjustable backrest, epoxy-coated steel frame and easy-clean vinyl upholstery. Sturdy, comfortable patient couch for clinics, GP rooms and triage areas.

The Vital Medical Ratchet-Type Examination Couch is a versatile, hard-wearing patient couch built for comfort, stability and everyday clinical use. Its backrest adjusts smoothly and locks securely via a 30 x 8mm ratchet mechanism with 25mm tubular support, allowing optimal patient positioning and easier caregiver access. The frame is fabricated from durable 32mm tubular mild steel finished in a white epoxy coating that resists corrosion and daily wear, while high-grade vinyl with a reinforced fabric backing covers all exposed surfaces for patient comfort and quick, hygienic cleaning. Four non-slip white rubber ferrules keep the couch firmly in place without marking or staining vinyl floors. Ideal for GP practices, hospital examination and admission rooms, triage and emergency units, wellness and physiotherapy centres, and first-aid stations.

Specifications

- Height: 820 mm
- Width: 610 mm
- Length: 1830 mm
- Frame: Epoxy-coated mild steel (32mm tubular)
- Backrest: Laser-cut ratchet mechanism with tubular support
- Upholstery: High-grade vinyl with fabric backing
- Feet: 4 ? white rubber ferrules (non-slip, non-marking)

SKU: VM-EXCRAT