



The medium stainless steel instrument trolley features a fully welded stainless steel frame and dual stainless steel shelves – upper and lower – mounted on four 75mm swivel castors. At 458 x 610 x 850mm, it offers a wider working surface than the small model, making it well suited to procedure rooms, theatres, and clinical environments where additional instrument capacity is needed without a large floor footprint.

This medium instrument trolley is fabricated in stainless steel with two 610 x 458mm shelves on swivel castors, giving a fully corrosion-resistant, easy-to-sanitise mobile surface. The seamless stainless construction withstands rigorous disinfection, making it well suited to theatres, treatment rooms and sterile areas that demand the highest hygiene standards.

SKU: VM-ITMSS