



The Vital Medical ILS Paramedic Bag (VM-ILSBC) is a fully stocked Intermediate Life Support bag providing advanced pre-hospital emergency care supplies. Adds IV cannulation (6 gauge sizes), larger syringes, and infant resuscitator over the BLS bag for ILS-level practitioner use.

What Is the ILS Paramedic Bag?

The VM-ILSBC (FAK-ILS) is a comprehensive Intermediate Life Support paramedic bag designed for ILS-level practitioners ? advanced emergency care technicians (AECTs) and emergency medical practitioners ? who are authorised to perform IV cannulation and administer fluids and medications at the scene and in transit. The ILS bag provides all BLS capabilities and adds the IV access, larger syringes, antiseptic solution, and infant resuscitator required for the expanded ILS scope of practice.

Contents and ILS-Specific Additions

The VM-ILSBC includes all standard first aid, wound care, airway, oxygen delivery, obstetric, and monitoring items found in the BLS bag, with the following ILS additions: full IV cannula range (14g, 16g, 18g, 20g, 22g, 24g ? 4 each), hypodermic needles 18g and 20g, syringes 10ml and 20ml (versus 5ml in BLS), antiseptic solution 50ml, rescue scissors 190mm (larger than BLS), and a

dedicated infant resuscitator in addition to the adult resuscitator. Complete wound care: bandages (50/75/100mm conforming, crepe, triangular), burn dressings 100mm ? 4, first aid dressings No.2/3/5, gauze, cotton wool, plasters, linen savers. Monitoring: pocket aneroid BP meter, stethoscope, pupil torch, digital thermometer. Obstetric: cord clamps ? 3, mucus extractor, maternity delivery pads.

Bag Design

The ILS bag features a distinctive red exterior with black straps, multiple internal compartments colour-coded for airway, vascular, monitoring, and obstetric items. Backpack-style carry straps and a rigid internal frame maintain bag shape and allow rapid one-handed access to individual compartments at the scene.

Technical Specifications

Specification	Detail
Specification	Value
Code	FAK-ILS
Level	Intermediate Life Support (ILS)
IV Access	Cannulas 14g?24g (4 each), needles 18g + 20g, syringes 10ml + 20ml
Airway	8? Guedel airways (000?5), 1? CPR mouthpiece
Oxygen Delivery	2? adult + 2? paediatric O2 masks + nebuliser masks
Resuscitation	1? adult + 1? infant resuscitator
Monitoring	Aneroid BP meter, stethoscope, pupil torch, thermometer
Obstetric	3? cord clamps, 1? mucus extractor, 2? maternity pads
Wound Care	Bandages, dressings, burn dressings, antiseptic 50ml
Application	ILS pre-hospital emergency, AECT scope of practice

SKU: VM-ILSBC