

The Vital Medical Professional Patella Hammer (VM-PHP) is a high-quality neurological reflex hammer with a durable rubber head and chrome-plated steel handle. Designed for accurate deep tendon reflex testing in clinical and hospital settings.

What Is the Professional Patella Hammer?

The VM-PHP is a classic Taylor-style reflex hammer built to professional quality standards for reliable daily use in clinical practice. The teardrop-shaped rubber head provides the precise percussive force required to elicit consistent deep tendon reflexes, while the balanced chrome-plated steel handle ensures comfortable grip during neurological examinations.

Clinical Uses and Benefits

The professional patella hammer is an essential diagnostic tool for assessing the patellar, brachioradialis, biceps, triceps, and Achilles reflexes. Graded reflex response (0?4+) provides important information about both upper and lower motor neuron integrity. This hammer's reliable rubber composition does not degrade with frequent clinical use, making it ideal for busy GP practices, emergency departments, and outpatient neurology clinics.

Build Quality and Maintenance

Constructed from chrome-plated steel with a quality vulcanised rubber head, the VM-PHP withstands the demands of daily clinical use. The rubber head resists hardening and cracking under normal storage conditions. The handle can be cleaned with standard surface disinfectants between patient encounters. Its classic design is familiar to all clinicians trained in standard neurological examination technique.

Technical Specifications

Specification	Detail
Specification	Value
Type	Taylor (teardrop) reflex hammer
Head Material	Vulcanised rubber
Handle Material	Chrome-plated steel
Style	Professional quality
Function	Deep tendon reflex assessment
Application	Neurology, GP, emergency
Sterilisation	Surface disinfectant wipe

SKU: VM-PHP