



Disposable peak flow mouthpieces supplied in packs of 100, available in white or blue. With an inner diameter of 24mm and outer diameter of 30mm, they are compatible with adult and paediatric peak flow meters and spirometers. CE and ISO approved for single-patient hygienic use in clinical and home settings.

Universal Fit for Adult and Paediatric Peak Flow Meters

These standard-dimension disposable mouthpieces are designed to fit adult and child peak flow meters ? including the Rossmax PF120A and PF120C ? as well as most spirometers using a compatible port size. With an inner diameter of 24mm and outer diameter of 30mm at a length of 5cm, they seat securely without requiring adapters, ensuring an airtight seal for accurate flow measurement. Available in white or blue, colour coding can be used to differentiate patient groups or devices in multi-patient environments.

Single-Use Hygiene for Multi-Patient Clinical Settings

Single-patient disposable use eliminates cross-contamination risk between patients ? a requirement for any shared peak flow meter used in GP rooms, asthma clinics, occupational health assessments, school health programmes, and hospital wards. Using a fresh mouthpiece per patient also removes the need for between-patient cleaning of the mouthpiece itself, streamlining workflow in high-volume environments without compromising infection control standards.

CE and ISO Approved ? Pack of 100

Each packet contains 100 mouthpieces (pack dimensions 410 ? 280 ? 90mm, 930g), providing sufficient stock for extended clinical use without frequent reordering. CE and ISO quality approvals confirm that the mouthpieces meet applicable medical device standards for material safety and dimensional consistency. Their straightforward design and low per-unit cost make them a practical consumable for both bulk clinical procurement and individual patient home monitoring kits.

Technical Specifications: