



The Silicone Adult Resuscitator CP (VM-RSSASCP) is the CP (Carry Pack) version of the Vital Medical silicone adult resuscitator, supplied with masks and accessories in a dedicated carry bag for organised storage and rapid deployment. It has a bag volume of 1500ml (± 200 ml), a stroke volume of ≥ 600 ml, a pressure limiting valve (40-60cmH₂O) and a 2000ml oxygen reservoir delivering an

FiO₂ of $\geq 85\%$ at a 15L/min oxygen flow rate. The ISO5356-1 22/15mm patient connector is compatible with adult masks, ET tube connectors and supraglottic airway devices. It measures 212 × 131mm and is autoclavable for multi-patient reuse. Available in South Africa from Vital Medical.

What Is the Silicone Adult Resuscitator CP?

The VM-RSSASCP is the CP (Carry Pack) version of the Vital Medical silicone adult resuscitator, supplied with masks and accessories in a dedicated carry bag for organised storage and rapid deployment at the point of need. The CP series meets the clinical performance requirements of the standard silicone resuscitator but is designed for paramedic kits, emergency trolleys and departments where a complete, ready-to-use resuscitator kit is required to be immediately available without individual component assembly.

Technical Performance

The adult CP resuscitator delivers a bag volume of 1500ml (± 200 ml) with a stroke volume of ≥ 600 ml, sufficient for adult tidal volumes in CPR (recommended 600ml at 12 BPM). The pressure limiting valve (40-60cmH₂O) prevents barotrauma, and the 2000ml oxygen reservoir ensures FiO₂ $\geq 85\%$ at a supplementary oxygen flow rate of 15L/min. The ISO5356-1 22/15mm patient connector is universally compatible with adult masks, ET tube connectors and supraglottic airway devices.

Infection Control, Reusability and Availability in South Africa

Silicone construction allows full autoclave sterilisation between patients, making the VM-RSSASCP a cost-effective choice for high-usage resuscitation environments. The carry bag maintains component organisation and protects the equipment from environmental contamination during storage. All components should be inspected before each use for cracks, valve patency and bag integrity. Available in South Africa from Vital Medical, the VM-RSSASCP suits paramedic kits, emergency trolleys and departments needing an immediately deployable adult resuscitator. Facilities can request a free quote through the Vital Medical website and receive a personalised response within 24 hours.

TECHNICAL SPECIFICATIONS TABLE

Specification	Detail
Brand	Vital Medical
SKU	VM-RSSASCP
Type	Manual BVM resuscitator – silicone – adult CP (in carry bag)
Patient Size	Adult (>40kg)
Bag Volume	1500ml ± 200ml
Stroke Volume	≥600ml
Reservoir Volume	2000ml
Dimensions	212 × 131mm
Pressure Limiting Valve	40-60cmH ₂ O
Patient Connector	ISO5356-1: Ø22/15mm
Sterilisation	Autoclavable
Supplied	In carry bag with masks

SKU: VM-RSSASCP