



The Charder Hebe 1 Digital Floor Scale (VM-SCH1), model SCAH1, is a slim, lightweight digital bathroom and home health scale with built-in BMI calculation for regular personal weight and health monitoring. It has a capacity of 150kg with 100g graduation, a 310 × 320mm platform, and an overall height of 25mm at 1.1kg. The BMI function stores up to four individual user height settings

and calculates body mass index from stored height and measured weight. It has a 1.0" LCD display, automatic power-off after 6 seconds, kg/lb interchangeable units, and runs on two 3V lithium batteries. It is CE certified. Available in South Africa from Vital Medical.

Slim Digital Floor Scale with BMI Function for Home Health Monitoring

The Charder Hebe 1 (VM-SCH1, SCAH1) is a slim, lightweight digital bathroom and home health scale with built-in BMI calculation, designed for regular personal weight and health monitoring in domestic and light clinical settings. The BMI function calculates body mass index from the user's stored height and current measured weight, displaying the result alongside the body weight reading. This allows the user to track their weight in the context of their height, which is the standard clinical measure of body weight status used in health assessment and chronic disease management. The Hebe 1 stores up to four individual user height settings, making it practical for households with multiple users who each want their BMI displayed after weighing without manually re-entering height each time.

Slim Profile and Lightweight Design

At 25mm high and 1.1kg, the Hebe 1 is one of the slimmest and lightest scales in the Charder range, easy to store under a bathroom vanity or step on and off without a raised platform edge. The 310 × 320mm platform provides a stable standing surface for standard adult users. The 1.0" LCD display is clearly readable in standard bathroom lighting conditions. Automatic power-off after 6 seconds conserves battery life, and the scale is powered by two 3V lithium batteries, which have a longer shelf life and more stable discharge characteristics than alkaline batteries in domestic conditions.

South African Home Health and Clinic Scale

The Hebe 1 is CE certified and suitable for South African home users, corporate wellness programmes and clinic waiting room weight monitoring applications requiring a lightweight, compact digital scale with BMI calculation for routine weight and health status monitoring. Available in South Africa from Vital Medical, the Hebe 1 is supplied to South African pharmacies, health clinics and home users seeking a reliable, slim digital floor scale with BMI function. Facilities can request a free

quote through the Vital Medical website and receive a personalised response within 24 hours.

TECHNICAL SPECIFICATIONS TABLE

Specification	Detail
Brand	Charder
SKU	VM-SCH1
Model	Hebe 1 (CODE: SCAH1)
Type	Digital floor scale with BMI
Capacity	150kg, 100g graduation
Overall Dimensions	310 × 320 × 25mm
Display	1.0" LCD
Power Supply	3V × 2 lithium batteries
Unit Weight	1.1kg
BMI Function	Yes (stores 4 user height settings)
Auto Power Off	After 6 seconds
Units	kg/lb interchangeable
Quality Standard	CE

SKU: VM-SCH1