



The Charder MS4300 Digital Baby Hanging Scale with Sling (VM-SBHMS4300) is a digital hanging baby scale designed for rapid, hygienic infant weight measurement using a suspended fabric sling, hung from any appropriate overhead point. It has a 25kg capacity with 10g resolution, a full-range tare function, a compact 9.5 × 16cm scale body and a 45 × 31cm sling, a 0.5" LCD display showing KG/LB/OZ, and a stabilisation time of 1-2 seconds. It weighs 265g gross, operates from 5°C to 35°C, has automatic external calibration, and runs on a 9V alkaline battery. It is CE certified and supplied complete with the scale unit, hook and baby sling. Available in South Africa from Vital Medical.

SHORT DESCRIPTION

The Charder MS4300 Digital Baby Hanging Scale with Sling (VM-SBHMS4300) is a digital hanging baby scale designed for rapid, hygienic infant weight measurement using a suspended fabric sling, hung from any appropriate overhead point. It has a 25kg capacity with 10g resolution, a full-range tare function, a compact 9.5 × 16cm scale body and a 45 × 31cm sling, a 0.5" LCD display showing KG/LB/OZ, and a stabilisation time of 1-2 seconds. It weighs 265g gross, operates from 5°C to 35°C, has automatic external calibration, and runs on a 9V alkaline battery. It is CE certified and supplied complete with the scale unit, hook and baby sling. Available in South Africa from Vital Medical.

LONG DESCRIPTION

Digital Baby Hanging Scale with Sling for Rapid Infant Weighing

The Charder MS4300 (VM-SBHMS4300) is a digital hanging baby scale designed for rapid, hygienic infant weight measurement using a suspended sling, a weighing method widely used in community health, nutrition and humanitarian settings where a flat-surface baby scale is impractical. The hanging scale with sling method allows the infant to be placed in a fabric sling that suspends from the hook of the scale, removing the need for a flat platform and allowing the scale to be hung from any appropriate overhead point, including door frames, tree branches, roof beams or purpose-made hanging stands in temporary clinical and community settings.

Community Health and Outreach Infant Weighing

The 25kg capacity and 10g resolution provide appropriate precision for infant weight monitoring from newborn through toddler age, and the full-range tare function allows the sling weight to be zeroed out so that only the infant's weight is displayed. The compact 9.5 × 16cm scale body and 265g gross weight make the MS4300 one of the most portable infant weighing instruments available. It can be carried in a bag alongside other clinical outreach supplies and deployed in any setting where a hanging point is available. The 9V alkaline battery provides reliable power in settings without mains electricity.

South African Community Health and Nutrition Programme Scale

The MS4300 is CE certified and is supplied complete with the scale unit, hook and baby sling, everything needed for immediate use. Available in South Africa from Vital Medical, the MS4300 is supplied to South African community health clinics, nutrition screening programmes, paediatric outreach teams, NGOs and child health monitoring programmes requiring a portable, easy-to-use infant hanging scale for community-based weight monitoring. Facilities can request a free quote through the Vital Medical website and receive a personalised response within 24 hours.

TECHNICAL SPECIFICATIONS TABLE

Specification	Detail
Brand	Charder
SKU	VM-SBHMS4300
Model	MS4300
Type	Digital baby hanging scale with sling
Capacity	25kg, 10g resolution
Scale Dimensions	9.5 × 16cm
Sling Dimensions	45 × 31cm
Display	0.5" LCD
Tare Range	Full range
Power Supply	9V alkaline battery
Units	KG/LB/OZ
Stabilisation Time	1-2 seconds
Operating Temperature	5°C to 35°C
Calibration	Automatic external
Gross Weight	0.265kg
Accessories	Scale × 1, Hook × 1, Baby Sling × 1
Quality Standard	CE

SKU: VM-SBHMS4300