



The Stethoscope Child Professional Dual Head (VM-STHCN1) is the Hi-Care Classic Dual Head Professional Child Stethoscope (STHPHC), a paediatric stethoscope with a stainless steel dual head, 37mm child head diameter, stainless steel retaining ring, binaural tubing and non-chill rings. The tube length is approximately 52cm and the complete stethoscope weighs 250g. The smaller head

gives good contact with the chest wall of infants and children, and the non-chill rings keep the chest piece comfortable against the skin during examination. It is CE and ISO approved. Available in South Africa from Vital Medical.

What Is a Paediatric Dual Head Stethoscope?

A paediatric stethoscope is an auscultation instrument sized for the smaller chest wall of infants and children. The Hi-Care Classic STPHC (VM-STHCN1) has a 37mm head diameter, smaller than the standard 45mm adult head. This gives consistent contact with the paediatric chest wall from infancy through school age. It also avoids the acoustic leakage that occurs when an adult-size head overhangs a small child's chest. This matters clinically because children's breath and heart sounds are higher in frequency and softer than adults', so accurate auscultation depends on good contact between the chest piece and the skin.

Construction, Comfort and Handling

The STPHC has a stainless steel dual head and stainless steel retaining ring, chosen for durability and easy cleaning in daily clinical use. Binaural tubing of approximately 52cm connects the head to the earpieces. Non-chill rings are fitted to the chest piece, which helps in paediatric examination because a cold chest piece can distress a child and cause movement that interferes with listening. At 250g, the stethoscope is lightweight enough for extended paediatric ward rounds. It carries CE and ISO approvals.

Paediatric Stethoscope Supply for South African Healthcare

The STPHC child professional stethoscope suits South African paediatric wards, general practitioners, community health clinics and nurses who regularly examine infant and child patients. Available in South Africa from Vital Medical, the VM-STHCN1 can be requested as a free quote through the Vital Medical website for healthcare facilities, paediatric departments and individual clinicians who need a correctly sized, CE- and ISO-approved professional stethoscope for daily use.

TECHNICAL SPECIFICATIONS TABLE

Specification	Detail
Brand	Hi-Care
Model	STHPHC (Hi-Care Classic Dual Head Professional)
SKU	VM-STHCN1
Type	Dual head stethoscope, child/paediatric
Head Material	Stainless steel
Head Diameter (Child)	37mm
Tube Length	Approx. 52cm
Retaining Ring	Stainless steel
Tubing	Binaural tubing assembly
Non-Chill Rings	Yes
Weight	250g
Quality Approvals	CE, ISO

SKU: VM-STHCN1