



The Professional Dual Head Adult Stethoscope (Classic Type) is a reliable, clinical-grade acoustic stethoscope designed for everyday diagnostic use. Featuring a dual-head chestpiece with both diaphragm and bell, it allows clinicians to assess a broad range of body sounds ? from high-frequency lung and bowel sounds to low-frequency heart murmurs ? making it a practical essential

for any healthcare professional.

The **Professional Dual Head Adult Stethoscope ? Classic Type** is a dependable, clinical-quality acoustic stethoscope built for routine diagnostic use across a wide range of healthcare settings. Designed for adult patients, it delivers the acoustic performance and durability that doctors, nurses, paramedics, and allied health professionals rely on daily.

Dual Head Chestpiece for Comprehensive Auscultation

The dual-head chestpiece is the defining feature of this stethoscope, providing two distinct auscultation surfaces in a single instrument:

- **Diaphragm** ? the larger, flat surface used to detect high-frequency sounds including breath sounds, normal heart sounds, and bowel sounds
- **Bell** ? the smaller, concave surface used to detect low-frequency sounds such as heart murmurs, extra heart sounds (S3, S4), and vascular bruits

This dual capability makes the stethoscope suitable for a broad range of clinical assessments without requiring the clinician to switch instruments.

Built for Clinical Reliability

The classic type design prioritises straightforward, consistent acoustic performance with a robust build suited to the demands of daily clinical use. Whether used in a busy hospital ward, a GP consulting room, an emergency setting, or during patient home visits, the Professional Dual Head Stethoscope delivers the clear, reliable sound transmission clinicians need to make confident

diagnostic decisions.

Suitable For:

General practitioners, nurses, paramedics, physiotherapists, medical students, hospital wards, emergency departments, outpatient clinics, and home visits.

Technical Specifications

Specification	Detail
Type	Professional acoustic stethoscope
Chestpiece	Dual head ? diaphragm and bell
Patient Size	Adult
Model	Classic Type

SKU: VM-STCC1