



The Surgical Set Basic 24pc (VM-SSB24) is a 24-piece stainless steel minor surgery instrument set supplied with a stainless steel instrument box (25 × 12 × 6cm). It contains 4 Bakhaus towel clamps, 2 Chaput tissue forceps, 6 Halsted mosquito forceps, Pean artery forceps, dressing forceps, Collin tissue forceps, standard tissue forceps, a 100ml gallipot, a Mayo-Hegar needle holder, a double-ended probe, a pair of Farabeuf baby retractors, Metzenbaum and Mayo curved scissors, and a #4 scalpel handle. All instruments are surgical-grade stainless steel and autoclavable at 134°C prevacuum or 121°C gravity. The set suits incision and drainage, excision of skin lesions, wound exploration and minor soft tissue surgery in GP procedure rooms, outpatient clinics and district hospitals. Available in South Africa from Vital Medical.

What Is the Basic Surgery Set 24pc?

The VM-SSB24 is a 24-piece general purpose minor surgical instrument set. It is designed for procedures such as incision and drainage, excision of skin lesions, wound exploration and minor soft tissue surgery. The instruments are pre-selected and supplied together in a stainless steel instrument box, so facilities do not need to source each instrument separately. This gives a consistent, standardised setup for minor surgical procedures and simplifies procurement and stock control.

Instrument Contents and Clinical Use

The set contains a stainless steel instrument box (25 × 12 × 6cm), four Bakhaus towel clamps (10.5cm) and two Chaput tissue forceps (14cm, 2×3 teeth). It also has six Halsted mosquito forceps (12.5cm, curved), a Pean artery forceps (14cm, straight), standard dressing forceps (14.5cm, straight), Collin tissue forceps (16cm, oval jaws) and standard tissue forceps (14.5cm, 1×2 teeth, straight). Suturing and dissection are covered by a Mayo-Hegar needle holder (18cm), a Metzenbaum curved scissor (14cm), a Mayo curved scissor (14cm) and a #4 scalpel handle. A double-ended probe (14.5cm), a pair of Farabeuf baby retractors (12cm) and a 100ml stainless steel gallipot complete the set. Together these instruments cover clamping, grasping, cutting, retracting, probing and holding fluids during a procedure.

Sterilisation, Setting and Procurement

All instruments in the VM-SSB24 are surgical-grade stainless steel and autoclavable at 134°C prevacuum or 121°C gravity, so the set can be reprocessed and used repeatedly. The instrument box serves as both storage and sterile presentation tray for the procedure. The set is suitable for GP procedure rooms, outpatient surgical clinics, district hospitals and private surgical suites performing minor to intermediate surgical procedures under local anaesthesia. Available in South Africa from Vital Medical. Healthcare facilities can request a free quote through the Vital Medical website and receive a personalised response within 24 hours.

TECHNICAL SPECIFICATIONS TABLE

Specification	Detail
Brand	Vital Medical
SKU	VM-SSB24
Set Contents	24 pieces
Instrument Box	Stainless steel, 25 × 12 × 6cm
Bakhaus Towel Clamps	4 × 10.5cm
Chaput Tissue Forceps	2 × 14cm (2×3 teeth)
Halsted Mosquito Forceps	6 × 12.5cm, curved
Pean Artery Forceps	14cm, straight
Dressing Forceps	Standard, 14.5cm, straight
Collin Tissue Forceps	16cm, oval jaws
Tissue Forceps	Standard, 14.5cm (1×2 teeth), straight
Gallipot	100ml, stainless steel
Needle Holder	Mayo-Hegar, 18cm
Probe	Double-ended, 14.5cm
Retractors	Farabeuf baby retractor, 12cm (pair)
Scissors	Metzenbaum curved 14cm; Mayo curved 14cm
Scalpel Handle	#4
Material	Surgical-grade stainless steel
Sterilisation	Autoclavable: 134°C prevacuum or 121°C gravity
Application	Minor/intermediate surgery, wound procedures
Setting	GP rooms, outpatient clinics, district hospitals, private surgical suites

SKU: VM-SSB24