



The Large Stainless Steel Dressing Trolley is a hand-polished, clinical-grade mobile trolley built from 25mm round tube stainless steel, featuring two stainless steel shelves, a 10mm chrome-plated guard rail on three sides, and four 75mm swivel castors. At 915mm wide, it provides an expansive, hygienic workspace for dressing procedures, instrument organisation, and sterile field preparation in

theatres, procedure rooms, and clinical wards.

The **Large Stainless Steel Dressing Trolley** is a premium-grade, fully mobile clinical trolley designed for use in operating theatres, procedure rooms, sterile environments, and hospital wards requiring a wide, hygienic working surface for dressing changes, minor surgical procedures, and instrument arrangement. Manufactured from 25mm round tube stainless steel with a hand-polished finish, it combines superior corrosion resistance with the clean, professional aesthetic required in sterile clinical environments.

Two Stainless Steel Shelves for Organised Clinical Use

The trolley features two stainless steel shelves – one at the top working level and one at a lower tier – providing organised, tiered storage and preparation space for dressings, instruments, sterile supplies, and procedure consumables. Both shelves are accommodated by four 16mm round tubes welded at the top and 220mm from the bottom of the frame respectively, creating a stable, secure shelf support system.

Chrome-Plated Guard Rail for Added Safety

A 10mm chrome-plated guard rail is fitted around three sides of the top frame, preventing instruments, dressings, and equipment from sliding off the working surface during procedures – particularly important when the trolley is in motion or positioned on uneven flooring. The rail maintains a clean, open front for easy access during clinical use.

Smooth, Controlled Mobility

The trolley is mounted on four 75mm swivel castors, enabling smooth, effortless movement across ward and theatre floors. The swivel action allows full 360° manoeuvring, making it easy to position precisely at the point of care without disrupting the clinical environment.

Hand-Polished Stainless Steel Throughout

The entire frame and shelving are constructed from stainless steel ? non-porous, corrosion-resistant, and fully compatible with hospital-grade disinfectants and steam cleaning. The hand-polished finish resists bacterial adhesion and maintains a hygienic surface even under frequent and rigorous decontamination, meeting the standards required for use in sterile clinical environments.

Suitable For:

Operating theatres, procedure rooms, sterile environments, dressing rooms, hospital wards, day surgery units, outpatient clinics, and any clinical setting requiring a large, hygienic mobile working surface.

Technical Specifications

Specification	Detail
Size	Large
Frame Material	25mm round tube stainless steel
Finish	Hand-polished stainless steel
Shelves	2 x stainless steel
Shelf Support	16mm round tubes welded at top and 220mm from bottom
Guard Rail	10mm chrome-plated, 3 sides
Castors	4 x 75mm swivel
Height	850mm
Width	915mm
Depth	458mm
Available Sizes	Small (458 x 458 x 850mm) / Medium (458 x 610 x 850mm) / Large (458 x 915 x 850mm)

SKU: VM-DTLSS