



The Medium Stainless Steel Dressing Trolley is a hand-polished, clinical-grade mobile trolley built from 25mm round tube stainless steel, featuring two stainless steel shelves, a 10mm chrome-plated guard rail on three sides, and four 75mm swivel castors. At 610mm wide, it strikes the ideal balance between workspace and footprint ? suited to a wide range of clinical environments from ward dressing bays to procedure rooms.

The **Medium Stainless Steel Dressing Trolley** is a versatile, premium-grade mobile clinical trolley designed for hospitals, clinics, procedure rooms, and ward environments where a practical mid-sized working surface is required. Manufactured from 25mm round tube stainless steel with a hand-polished finish, it offers a wider working area than the small model while remaining compact enough for use in most standard clinical spaces ? making it the most adaptable size in the range.

Two Stainless Steel Shelves for Organised Clinical Use

The trolley features two stainless steel shelves ? one at the top working level and one at a lower tier ? providing organised, tiered storage and preparation space for dressings, instruments, sterile supplies, and procedure consumables. Both shelves are supported by four 16mm round tubes welded at the top and 220mm from the bottom of the frame, ensuring a stable, secure structure under clinical use.

Chrome-Plated Guard Rail for Added Safety

A 10mm chrome-plated guard rail is fitted around three sides of the top frame, preventing instruments, dressings, and equipment from sliding off the working surface during procedures or while the trolley is repositioned. This is an essential clinical safety feature across all sizes of the range.

Smooth, Controlled Mobility

Four 75mm swivel castors provide smooth, effortless movement across ward and theatre floors. The

610mm width gives the trolley a stable, well-balanced base that moves confidently across a range of clinical floor surfaces while remaining easy to position precisely at the bedside or procedure table.

Hand-Polished Stainless Steel Throughout

The fully stainless steel frame and shelving is non-porous, corrosion-resistant, and fully compatible with hospital-grade disinfectants and steam cleaning. The hand-polished finish resists bacterial adhesion and maintains a hygienic surface under the frequent and rigorous decontamination required in sterile clinical environments.

Suitable For:

Hospital wards, dressing rooms, procedure rooms, outpatient clinics, theatre ante-rooms, day surgery units, consulting rooms, and any clinical setting requiring a mid-sized, hygienic mobile dressing trolley.

Technical Specifications

Specification	Detail
Size	Medium
Frame Material	25mm round tube stainless steel
Finish	Hand-polished stainless steel
Shelves	2 x stainless steel
Shelf Support	16mm round tubes welded at top and 220mm from bottom
Guard Rail	10mm chrome-plated, 3 sides
Castors	4 x 75mm swivel
Height	850mm
Width	610mm
Depth	458mm
Available Sizes	Small (458 x 458 x 850mm) / Medium (458 x 610 x 850mm) / Large (458 x 915 x 850mm)



Trolley Dressing Medium Stainless Steel 610 x 458mm

SKU: VM-DTMSS