



The Small Stainless Steel Dressing Trolley is a compact, hand-polished clinical-grade mobile trolley built from 25mm round tube stainless steel, featuring two stainless steel shelves, a 10mm chrome-plated guard rail on three sides, and four 75mm swivel castors. At 458mm wide, it is ideal for smaller clinical spaces, treatment rooms, and consulting rooms where a hygienic, space-efficient dressing trolley is required.

The **Small Stainless Steel Dressing Trolley** is a compact, premium-grade mobile clinical trolley designed for use in consulting rooms, treatment areas, procedure rooms, and clinical wards where space is limited but hygiene and clinical functionality cannot be compromised. Manufactured from 25mm round tube stainless steel with a hand-polished finish, it delivers the same corrosion-resistant, sterile-environment performance as its larger counterparts in a footprint suited to smaller clinical spaces.

Two Stainless Steel Shelves for Organised Clinical Use

Despite its compact size, the trolley features two stainless steel shelves – one at the top working level and one at a lower tier providing organised, tiered storage and preparation space for dressings, instruments, sterile supplies, and procedure consumables. Both shelves are supported by four 16mm round tubes welded at the top and 220mm from the bottom of the frame, creating a stable, secure shelf structure.

Chrome-Plated Guard Rail for Added Safety

A 10mm chrome-plated guard rail fitted around three sides of the top frame prevents instruments and dressings from sliding off the working surface during procedures or while the trolley is in motion – an essential safety feature in busy clinical environments regardless of trolley size.

Smooth, Controlled Mobility

Four 75mm swivel castors provide smooth, effortless movement across ward and clinical room floors.

The compact 458mm x 458mm base makes the small dressing trolley particularly easy to manoeuvre in tighter spaces such as consulting rooms, treatment cubicles, and smaller procedure areas.

Hand-Polished Stainless Steel Throughout

The fully stainless steel frame and shelving is non-porous, corrosion-resistant, and compatible with hospital-grade disinfectants and steam cleaning. The hand-polished finish maintains a hygienic, easy-to-clean surface under frequent decontamination suitable for use in sterile clinical environments.

Suitable For:

GP consulting rooms, treatment rooms, procedure cubicles, nursing stations, outpatient clinics, day surgery units, and any clinical setting requiring a compact, hygienic mobile dressing trolley.

Technical Specifications

Specification	Detail
Size	Small
Frame Material	25mm round tube stainless steel
Finish	Hand-polished stainless steel
Shelves	2 x stainless steel
Shelf Support	16mm round tubes welded at top and 220mm from bottom
Guard Rail	10mm chrome-plated, 3 sides
Castors	4 x 75mm swivel
Height	850mm
Width	458mm
Depth	458mm
Available Sizes	Small (458 x 458 x 850mm) / Medium (458 x 610 x 850mm) / Large (458 x 915 x 850mm)

SKU: VM-DTSSS