



The FS903 Basic Steel Wheelchair (VM-WCBK903) is a steel-frame, PVC-upholstered adult wheelchair with an 18-inch (46cm) seat width, 100kg weight capacity, and detachable armrests and footrests. The folding X-frame collapses for storage and vehicle transport. Overall width 50cm, total height 89cm, packed dimensions 80×28×91cm, net weight 17.9kg. Supports both self-propulsion and attendant-assisted use

Durable Steel Frame Wheelchair for Hospital, Rehabilitation, and Home Care

The FS903 is a steel-frame, PVC-upholstered wheelchair designed for reliable daily use in hospital wards, rehabilitation centres, frail care facilities, and home care settings. The robust steel frame provides durable support under continuous use in demanding clinical environments, and the 100kg weight capacity covers the large majority of the adult user population. The 18-inch (46cm) seat width provides a comfortable seated position for standard adult body dimensions, and the PVC seat and back surfaces are durable and cleanable — a practical requirement in shared clinical and care environments. Rear push handles allow attendant-assisted propulsion, while the large rear wheels with handrims support self-propulsion for users with adequate upper body strength, making the FS903 suitable for both dependent and partially independent users.

Detachable Armrests, Footrests, and Folding Frame

The detachable armrests and footrests on the FS903 are a clinically significant design feature, allowing transfer to and from the wheelchair — to a bed, toilet, vehicle seat, or examination couch —

without requiring the user or carer to negotiate around fixed side components. Detachable armrests reduce lateral transfer distance and injury risk, and are particularly important for users with limited lower limb strength or those undergoing rehabilitation where transfer independence is a goal. The folding X-frame allows the chair to be collapsed quickly for storage in a compact space or loading into a standard passenger vehicle, with packed dimensions of 80×28×91cm and a net weight of 17.9kg. The 50cm overall width is compatible with standard doorway and hospital corridor clearances, making the FS903 practical for indoor use in both residential and clinical environments.

South African Supply for Healthcare, Rehabilitation, and Home Use

The FS903 is a cost-effective, dependable standard wheelchair suited to South African public and private hospitals, rehabilitation units, nursing homes, disability equipment loan programmes, and home care environments where a practical, functional adult wheelchair is required for daily mobility support. Available in South Africa from Vital Medical, the FS903 is supplied to healthcare facilities, NGOs, and individual users requiring a durable attendant or self-propelled steel wheelchair at a standard adult seat width.

TECHNICAL SPECIFICATIONS TABLE

Specification	Detail
Model	FS903
SKU	VM-WCBK903
Type	Standard Steel/PVC Wheelchair
Frame Material	Steel
Upholstery	PVC
Seat Width	18 inch (46cm)
Weight Capacity	100kg
Armrests	Detachable
Footrests	Detachable
Frame Style	Folding X-frame
Overall Width	50cm

Specification		Detail
Total Height	89cm	
Packed Dimensions	80 × 28 × 91cm	
Net Weight	17.9kg	
Gross Weight	20.9kg	

SKU: VM-WCBK903