



# Steel Furniture Catalogue 2022



### **MEDIUM DRESSING TROLLEY**

Stainless steel trolley, stainless steel shelves, 10mm chrome plated guard rail covering 3 sides, 4 x 75mm swivel castors



### **SMALL DRESSING TROLLEY**

Stainless steel trolley, stainless steel shelves, 10mm chrome plated guard rail covering 3 sides, 4 x 75mm swivel castors



### **LARGE DRESSING TROLLEY**

Stainless steel trolley, stainless steel shelves, 10mm chrome plated guard rail covering 3 sides, 4 x 75mm swivel castors



### **MEDIUM INSTRUMENT TROLLEY**

Epoxy coated trolley with stainless steel shelves, 4 x 25mm chrome plated caps to seal top tubes, 4 x 75mm swivel castors



### **LARGE INSTRUMENT TROLLEY**

Epoxy coated trolley with stainless steel shelves, 4x 25mm chrome plated caps to seal top tubes, 4 x 75mm swivel castors



### **SMALL INSTRUMENT TROLLEY**

Epoxy coated trolley with stainless steel shelves, 4 x 25mm chrome plated caps to seal top tubes, 4 x 75mm swivel castors



### **STAINLESS STEEL INSTRUMENT TROLLEY**

Stainless steel trolley with stainless steel shelves, 4 x 25mm chrome plated caps to seal top tubes, 4 x 75mm swivel castors



### **STAINLESS STEEL FIXED HEIGHT OVERBED TABLE**

L 1090 x W 400 x H 900mm, stainless steel frame and top, 4 x 75mm swivel castors



### **WOOD TOP FIXED HEIGHT OVERBED TABLE**

L 1090mm x W 400mm x H 900mm, mild steel epoxy coated frame, 4 x 75mm swivel castors, solid formica wood top



### **FIXED HEIGHT OVERBED TABLE**

OE 22

L 1090mm x W 400mm x H 900mm, epoxy powder coated mild steel frame and top, 4 x 75mm swivel castors



### **MAYO TABLE WITH LARGE TRAY**

Stainless steel, 4 x 50mm swivel castors, stainless steel tray, extension rod (800-1300mm), tray support



### **HI-LOW OVERBED TABLE – SPECIALISED STAINLESS STEEL TOP**

Epoxy powder coated mild steel base, specialised stainless steel top



### **HI-LOW OVERBED TABLE – SPECIALISED STAINLESS STEEL TOP**

Epoxy powder coated mild steel, formica wood top, lever handle on gas spring to adjust, 800mm-1100mm



### **HI-LOW OVERBED TABLE**

Epoxy powder coated mild steel, formica wood top, lever handle on gas spring to adjust, 800mm-1100mm



### **DOUBLE FOOT STOOL**

Mild steel epoxy powder coated, non-slip rubber top, aluminium edging, 4 x ferrules, 510mm x 660mm x 460mm



### **SINGLE FOOT STOOL**

Mild steel, epoxy powder coated, non-slip rubber top, aluminium edging, 4 x 25mm ferrules



### **DOUBLE FOOT STOOL**

Stainless steel with checkered aluminium top, 4x ferrules, 510mm x 660mm x 460mm



### **SINGLE FOOT STOOL**

Stainless steel, aluminium checkered top, 4 x 25mm ferrules



### **VISITORS BENCH**

Epoxy powder coated frame, available in 915mm, 610mm and 450mm, 4 legs fitted with ferrules, laminated wooden 20mm thick top, 30mm x 100mm centre hole, clear varnish finish



### **MISSION LOCKER**

Formica wood top, mild steel epoxy powder coated, locker with draw and 4x castors



### **COMBINATION BEDSIDE LOCKER – SHORT VERSION**

Epoxy powder coated, formica wood top, 1 x pull out drawer, 1 x cupboard with single roller catch, 1 x shelf inside cupboard, 1 x epoxy coated towel rail, combination locker



### **SINGLE BEDSIDE LOCKER**

Epoxy powder coated, formica wood top, 1 x pull out drawer, 1 x cupboard with single roller catch, 1 x shelf inside, 1 x epoxy coated towel rail



### **COMBINATION BEDSIDE LOCKER – TALL VERSION**

Epoxy powder coated, formica wood top, 1 x pull out drawer, 1 x cupboard with single roller catch, 1 x shelf inside cupboard, 1 x epoxy coated towel rail, 700mm x 430mm x 1060mm combination locker



### **SINGLE BEDSIDE LOCKER**

Epoxy powder coated, stainless steel top, 1 x pull out drawer, 1 x cupboard with single roller catch, 1 x shelf inside, 1 x epoxy coated towel rail



### **NEO-NATAL BED END TROLLEY**

Mild steel epoxy powder coated trolley, formica wood top, 1 drawer, 1 wire basket, 2 75ml brake castors, 2 plain castors



### **DOUBLE DOOR DRUG CABINET**

Epoxy powder coated, 1 fixed shelf, 1 adjustable shelf, lockable, schedule 7 lockable drug cabinet inside



### **LARGE POISON CABINET**

Epoxy powder coated, shelf inside, lockable door



### **LARGE DEADLY DANGEROUS CABINET**

Epoxy powder coated, shelf inside, lockable door, schedule 7 lockable drug cabinet inside



### **POISON CABINET**

Epoxy powder coated, shelf inside, lockable door



### **ONE ON TOP OF ANOTHER POISON CABINET**

Epoxy powder coated, 2 lockable doors



### **ANAESTHETIC TROLLEY - 1 DRAWER**

1 pull-out lockable drawer,  
stainless steel shelf,  
10mm chrome plated  
guard rail, 4 x 75mm  
swivel castors



### **ICU BED END TROLLEY**

Mild steel epoxy  
coated, formica wood  
top with aluminium  
edging, pull out top  
drawer - optional, 4  
drawers on left hand  
side 1 draw and 1 door  
on right hand side,  
needle box holder, soap  
dispenser, x ray holder



### **SPECIALISED 2 DRAWER ICU BED END TROLLEY**

Custom designed  
and manufactured  
to client  
specifications



### **5 DRAWER EMERGENCY TROLLEY**

Mild steel epoxy  
coated, stainless steel  
top, stainless steel top  
rail, 5 lockable drawers,  
drip stand holder,  
oxygen holder, CPR  
board, 2 x 125mm  
break castors 2 x  
125mm plain castors



### **ANAESTHETIC TROLLEY - 6 DRAWER 3 SHELVES**

6 drawers, 3 shelves,  
epoxy powder coated,  
stainless steel worktop  
surface, metal clips on  
either side, 4 x  
125mm castors



### **ANAESTHETIC TROLLEY - 6 DRAWER ONLY**

6 drawers, epoxy  
powder coated,  
stainless steel  
worktop surface, 4  
x 125mm castors,  
metal clips on  
either side



### **STAINLESS STEEL MOBILE DRIP STAND**

Stainless steel, 5 x 50mm swivel castors, adjustable T-rod, adjustable height (1300mm-2000mm)



### **CUSTOM MADE MOBILE DRIP STAND**

Epoxy powder coated, 5 x 50mm swivel castors, adjustable T-rod, adjustable height (1300mm-2000mm)



### **MOBILE DRIP STAND – 3 LEG**

Epoxy powder coated, 3 x 75mm swivel castors, adjustable T-rod, adjustable height (1300mm-2000mm)



### **MOBILE DRIP STAND – 5 LEG**

Epoxy powder coated, 5 x 50mm swivel castors, adjustable T-rod, adjustable height (1300mm-2000mm)



### **CUSTOMISED HEAVY BASE DRIP STAND**

Designed and manufactured according to client specifications



### **SPECIALISED MOBILE STAINLESS STEEL DRIP STAND**

Stainless steel brush finished, oxygen holder, adjustable height, 5 legged base, adjustable T-rod





### **OXYGEN CYLINDER STAND - UPRIGHT**

Size H 680 x L400  
x W400mm, epoxy  
coated mild steel,  
transports a  
180mm diameter  
cylinder



### **KICK-A-BOUT BUCKET TYPE**

Epoxy coated mild  
steel, single ring to  
hold bucket,  
3x50mm swivel  
castors



### **KICK-A-BOUT BUCKET TYPE**

Stainless steel,  
single ring to hold  
bowl, 3x50mm  
swivel castors



### **KICK-A-BOUT BOWL TYPE**

Stainless steel,  
single ring to hold  
bowl, 3 x 50mm  
swivel castors



### **DOUBLE BOWL STAND**

Stainless steel, H  
835mm x L  
610mm x W 390  
mm, 295mm  
diameter rings to  
hold bowls. Bowls  
sold separately



### **SINGLE BOWL STAND**

H 835mm x L  
390mm x W 390  
mm, stainless  
steel, single ring  
295mm diameter,  
4 swivel castors.  
Bowl sold  
separately



### **SOLID LINEN RECEIVER**

Epoxy powder coated, 610mm diameter ring, 610x900mm, canvas bag with nylon draw cord sold separately



### **26L CLINIBIN**

Epoxy powder coated, removable plastic bin with handle, bin lid activated by pedal, H 890mm x W 490mm x L 520mm



### **85L CLINIBIN**

Epoxy powder coated, removable plastic bin with handle, bin lid activated by pedal, H 1080mm x W 530mm x L 520mm



### **TWO DRAWER DRESSING TROLLEY**

Epoxy coated mild steel trolley and drawers, stainless steel shelf, 10mm chrome plated guard rail, 2 x pullout drawers



### **CUSTOMISED THREE TIER TROLLEY**

Epoxy coated silver vein, 3 tiers, 3 stainless steel trays, 940mm x 900mm x 530mm, specialised castors



### **SOLID LINEN RECEIVER WITH 3 SHELVES**

Epoxy powder coated, 445mm diameter ring, epoxy coated shelves 295mm apart, canvas bag with nylon draw cord sold separately



### **THREE TIER FOOD OR LINEN TROLLEY**

Mild steel epoxy coated, 3 tiers, 3 mild steel epoxy coated trays, 4x75mm swivel castors, 940mm x 900mm x 530mm



### **UTILITY DRUG CABINET**

Epoxy powder coated, 2 lockable doors, writing shelf, stainless steel tray with bowl holders, stainless steel tray with galvanized rail



### **PATIENT TROLLEY WITH STAINLESS STEEL TOP**

L 2000mm x 660mm x 660mm, epoxy powder coated frame, stainless steel top and stainless steel handles on each end, drip rod, oxygen cylinder holder, document/blanket tray, 2 break 2 swivel castors, 3 restraining straps



### **TWO SECTION HI-LO BED**

Two section hi-lo bed, manual gas lever backrest, hydraulic hi-lo foot pedal mechanism, removable head and foot boards, 125ml brake castors and 125ml normal castors



### **EXAMINATION COUCH WITH BREATHING HOLE**

Epoxy powder coated frame, ratchet operated backrest, vinyl covered chip foam (available in any colour), 120 x 120mm breathing hole



### **EXAMINATION COUCH**

Epoxy powder coated tube frame, ratchet operated backrest, vinyl covered chip foam (available in any colour)



### **FOLDABLE COUCH**

Specialised item,  
made on request



### **HOSTEL BED**

Specialised hostel beds  
made to order



### **CHILD COT**

Epoxy powder coated,  
sides can move up  
and down on movable  
guide, 2x brake  
castors, 2x swivel  
castors, 50x50x5mm  
mesh wire mattress  
base, mattress sold  
separately



### **SPECIALISED BASINET**

Specialised design,  
made to order.



### **CHILD COT**

Epoxy powder coated,  
sides can move up and  
down on movable guide,  
2x brake castors, 2x  
swivel castors,  
50x50x5mm mesh wire  
mattress base, mattress  
sold separately



### **BED SCREEN**

Epoxy powder  
coated, each panel  
is 1700mm x  
600mm, available  
in 3,4 & 5 fold.  
Curtain sold  
separately



### **SPECIALISED TEA TROLLEY**

Specialised item,  
made on request



### **TWO TIER TROLLEY**

Mild steel epoxy  
coated frame, 2  
tiers, 2 stainless  
steel trays, 4 x  
75mm swivel  
castors



### **BASINETTE**

L 700mm x W 440mm x H  
820mm, clear acrylic crib,  
4 x 50mm swivel castors,  
water resistant taffeta  
cover mattress



### **ANAESTHETIC STOOL**

Adjustable  
aluminium seat  
top, mild steel  
epoxy powder  
coated, 4 x 50mm  
swivel castors



### **OPERATION STOOL**

Epoxy coated mild  
steel, adjustable  
upholstered seat  
top, 4 x 50mm  
swivel castors



### **OPERATION STOOL**

Epoxy coated mild  
steel, adjustable  
upholstered seat  
top, 4 x ferrules



### **PATIEND CHART HOLDER**

Epoxy powder coated,  
2 pockets (width  
60mm each) 10mm  
round bar metal  
hooks covered with  
plastic protective  
tubing



### **STAINLESS STEEL PATIENT CHART HOLDER**

Stainless steel, 2  
pockets (width 60mm  
each) 10mm round  
bar metal hooks  
covered with plastic  
protective tubing.  
Made on request.



### **X RAY VIEWING BOX**

Epoxy powder  
coated light box to  
hold x-rays to view




### **LIFTING POLE (MONKEY POLE) WITH HANDLE**

Free standing,  
epoxy powder  
coated frame,  
triangular shaped  
handle with chain  
link, 4 holes on top  
horizontal bar to  
locate handle



### **URINE AND CATHETER BAG HOLDER**

Stainless steel  
urine bag holder



### **URINE AND CATHETER BAG HOLDER**

Mild steel, epoxy  
powder coated  
urine bag holder



### **OXYGEN CYLINDER CARRIER – SABS LEANING TYPE**

Epoxy coated mild steel, transports 180mm diameter cylinder, 2 front castors (150mm), 2 back castors (75mm), 10mm base plate



### **MOBILE MONITOR STANDS**

Specialised item, made on request



### **BED PAN AND URINAL RACK**

Stainless steel round bar made to hold bed pans and urinals, with stainless steel drip tray



### **INSTRUMENT CABINET**

Epoxy powder coated, lockable, double door, 2 to 5 adjustable glass shelves



### **INSTRUMENT CABINET**

Epoxy powder coated, lockable, double door, 2 adjustable glass shelves



**[www.vitalmedical.co.za](http://www.vitalmedical.co.za)**

**E: [info@vitalmedical.co.za](mailto:info@vitalmedical.co.za)**

**T: 087 265 8944**

**54 Augusta Road  
Johannesburg**