

The VM-EXG3 is a three-section gynaecological examination couch designed for pelvic, obstetric, and gynaecological examinations across GP consulting rooms, women's health clinics, antenatal care settings, and colposcopy suites. It features an independently adjustable backrest, an independently adjustable leg section for lithotomy positioning, and is supplied complete with stirrups, poles, and straps as standard. The couch measures 1830 x 610 mm with a fixed height of 760 mm, and is upholstered in a smooth, easy-clean clinical fabric suited to high-turnover gynaecological environments.

Three-Section Configuration for Precise Gynaecological Positioning

The VM-EXG3 three-section examination couch provides independent adjustment of both the backrest and leg sections, allowing clinicians to position patients precisely for the full range of gynaecological and obstetric procedures. The adjustable backrest supports patients in a semi-reclined or upright position for comfort during consultation, while the independently adjustable leg section enables the lithotomy position required for cervical smears, pelvic examinations, colposcopy, and minor gynaecological interventions. The ability to adjust both sections independently, rather than as a linked mechanism, gives clinicians the flexibility to optimise patient position for each procedure without compromising clinical access or patient comfort.

Stirrups, Poles, and Straps Included as Standard

The VM-EXG3 is supplied complete with stirrups, poles, and straps, providing the limb support required for lithotomy positioning without supplementary procurement. The inclusion of a full stirrup assembly makes this couch a ready-to-use solution for gynaecological examination rooms from the point of delivery, suitable for routine cervical screening, antenatal pelvic assessment, and family planning consultations. The upholstery is a smooth, easy-clean clinical fabric that can be wiped down quickly with standard hospital-grade disinfectants between patients, supporting infection control in high-throughput women's health settings.

Dimensions and Availability in South Africa

The VM-EXG3 has an overall length of 1830 mm, a width of 610 mm, and a fixed height of 760 mm, providing a stable and appropriately sized platform for adult gynaecological examination across a wide range of patient sizes. Available in South Africa from Vital Medical, the VM-EXG3 is suited to gynaecology and obstetric practices, GP consulting rooms, women's health clinics, family planning clinics, antenatal care settings, and colposcopy suites requiring a complete, purpose-built examination solution.

Technical Specifications

Specification	Value
Model	VM-EXG3
Couch Type	Three-section gynaecological examination couch
Sections	Backrest, seat, leg (3 sections)
Backrest	Independently adjustable
Leg Section	Independently adjustable
Upholstery	Easy-clean clinical fabric
Stirrups	Included
Poles	Included
Straps	Included
Length	1830 mm
Width	610 mm
Height	760 mm
Suitable Environments	Gynaecology, obstetrics, women's health, GP rooms, colposcopy, family planning, antenatal care

SKU: VM-EXG3